

LARYNGOSCOPE SET *SPECIAL OPERATIONS*

REV042420

The lightweight, compact, and highly innovative Special Operations Laryngoscope Set from North American Rescue is designed to operate effectively under the most extreme environmental conditions. Professionals rely on it. To ensure the laryngoscopes perform as long and as well as necessary, the sets are packaged with spare LED bulbs and batteries.

The fiber optic handle is water-proof, sand-proof, and dust-proof, and features a textured gripping surface. For convenience, the set includes single-use fiber optic blades made from an enhanced polymer for exceptional strength and rigidity. For added convenience, the handle's removable light source is designed to work with any standard (green-spec) fiber optic blade.

The complete set rolls up for easy storage and easy transport.

Kit Contents:

- ENT Pocket Light™ with two AAA Batteries
- Laryngoscope Handle
- Miller #3 Laryngoscope Blade
- Laryngoscope Pouch
- Spare LED Exam light bulb
- Two spare AAA batteries



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**

LARYNGOSCOPE SET *SPECIAL OPERATIONS*

REV042420

*Lightweight, compact
and highly innovative*

Features:

- Lightweight & compact kit rolls up for easy storage
- Resistant to dust, sand and water
- Operates in severe environmental conditions
- Textured, water-proof & dust-proof handle
- Removable 80 Lumen LED light source works with any standard fiber optic blade
- Single-use fiber optic blades of enhanced polymer for exceptional strength and rigidity

Specifications:

- Open: L 8 in. x W 5.5 in. x D 2 in.
- Closed: L 8 in. x W 2 in. x D 2 in.
- Weight: 10 oz

Available
Color:

BLK

Color	Item#	NSN#
BLK	10-0001	6515-01-450-9790



*Perform direct visual laryngoscopy
even in the most severe environments*



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH A

MISSION®