

IPOK[®]

Individual Patrol Officer Kit



Compact hemorrhage control kit designed to fit into a BDU pocket...

The Individual Patrol Officer Kit (IPOK) from North American Rescue is designed to provide personnel with a compact and durable individual hemorrhage control kit to treat bleeding from penetrating and other traumatic injuries.

These kits are packaged for small cube space and designed to fit into a BDU pocket, vest pouch, or individual bag, which allows personnel to keep a compact bleeding control kit on their person, where it is needed most. The contents are vacuum sealed in a rugged and durable packaging that is easy to open.

Features:

- Compact, vacuum-sealed hemorrhage control kit designed to be carried on your person where it is needed most
- Easy-to-open streamlined kit
- Designed to fit into a BDU pocket, vest pouch, or individual bag
- Your Choice of S-Rolled Gauze or Hemostatic

Specifications:

- Packaged: H 4 in. x W 6.5 in. x D 2.75 in.
- Weight: 7.4 oz



IPOK[®] Individual Patrol Officer Kit

Hemostatic/Dressing	Item#	NSN#
S-Rolled Gauze	80-0167**	
Combat Gauze LE	80-0168	



Kit Contents:

- 1 x pair Black Responders Gloves, lg.
- 1 x C-A-T[®], Black (Combat Application Tourniquet[®])
- 1 x Your Choice of Gauze (Wound Packing Gauze or Combat Gauze LE)
- 1 x Flat ETD[™] 4 in. Emergency Trauma Dressing



NORTH AMERICAN RESCUE[®]
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A MISSION[®]

IPOK® Dressing Options

REV101422



NAR Wound Packing Gauze *Option 1 for the IPOK® (Item #80-0167)*

- Soft, white gauze is 1-ply, 3 in. x 5 yd length
- Indicated for use as a wound packing gauze and includes embedded X-Ray detectable line
- Z-fold packaging configuration allows controlled application as a complete unit or in a continuous linear feed
- Used as a packing gauze for wounds, as a backing gauze for hemostatic agents, in conjunction with a compression bandage to control hemorrhage, or for minor wound bandaging



Combat Gauze LE Gauze *Option 3 for the IPOK® (Item #80-0168)*

- 3 in x 4 yd (12 ft)
- Z-Fold, Vacuum Packed for easy application
- Effective, Simple, and Safe – Does not generate heat
- Easy to use and remove
- Flexible and pliable

**Choose one
of these two
dressing options
to complete your
IPOK® Kit**



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®