

# NAR S-ROLLED GAUZE

REV052620

## **Compact, easy to use, S-Rolled cotton gauze designed for both linear and basic bandaging**

Made of premium cotton, S-Rolled Gauze from North American Rescue can be used to control hemorrhage in conjunction with a compression bandage, used as backing gauze for hemostatic dressings or for minor wound bandaging.

NAR S-Rolled Gauze is packaged in a unique dispenser that allows for controlled application and protects unused portions from contamination during application. Its low cube space requirement and rugged, durable vacuum sealed packaging make it easy to fit in both individual first aid and medic kits.

### **Features:**

- Unique and convenient packaging configuration and dispenser allows controlled application as a complete unit or in a continuous linear feed
- Unused portion remains protected from contamination during application
- Used as backing gauze for hemostatic agents, in conjunction with a compression bandage to control hemorrhage
- Vacuum-sealed packaging branded with NAR's signature Red Tip Technology® easy-open tear notches.

### **Specifications:**

- Vacuum Packed
- Packaged: L 3 in. x W 2 in. x D 1 in.
- Gauze: 6-ply, W 4.5 in. x L 4.1 yd
- Weight: 1.2 oz

### **NAR S-Rolled Gauze**

| Item#   | NSN#             |
|---------|------------------|
| 30-0003 | 6510-01-605-2234 |

CE  
0050

NOT MADE WITH  
NATURAL RUBBER LATEX

STERILE



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**