

QUIKCLOT™ Z-FOLD COMBAT GAUZE

REV060618



Stops Arterial and Venous bleeding in seconds...

Fast-acting Combat Gauze® is the ultimate wound care dressing. Tailored to the needs of combat and tactical medical personnel, this hemostatic dressing combines surgical gauze with a proprietary inorganic material to deliver three important benefits. Stops blood flow. Inert. Non-allergenic.

The material, pliable and familiar to medical personnel, can be fitted to any size or shape wound (including penetrating wounds). Plus, it is easily removed once clotting has taken place.

The Committee on Tactical Combat Casualty Care recently recommended Combat Gauze® as the first-line hemostatic product.

Features:

- Advanced hemostatic agent with the simplicity of gauze
- Easily removed once clotting occurs
- Can be fitted to any size or shape wound
- Non-allergenic, Latex-Free
- Flexible & pliable, able to contour to all wounds & does not generate heat
- Individually wrapped in easy-rip, military-grade foil pouch
- Featuring an embedded x-ray detectable strip



Standard Specifications:

- Gauze: L 12 ft. x W 3 in.
- Pkg: H 5 in. x W 5.5 in. x D 0.625 in.
- Weight: 0.8 oz.

XL Specifications:

- Gauze: 2-ply L 4 yd. x W 4 in.
- Pkg: H 6.5 in. x W 6.5 in. x D 0.75 in.
- Weight: 1.75 oz.

Size	Item#	NSN#
STD	30-0039	6510-01-562-3325
XL	30-0070	6510-01-600-4705



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH  **MISSION®**