

TACTICAL COMPRESSION WRAP



Most Versatile Bandage System Ever Developed:

When tactical medics pack their gear, they want minimal items with maximum application. They want the Tactical Compression Wrap from North American Rescue. A lightweight, compact wrap with contact bar and hook/loop attachment, this elastic bandage is easy to apply under the most stressful conditions.

An unsurpassed compression bandage, the Tactical Wrap can be used for extremity, head, shoulder and abdominal wounds, or for buddy splinting-and that's just the beginning of its circumferential applications.

Developed by a Special Forces Medic and has been tested with great success in both Afghanistan and Iraq. When space is at a premium, conditions are harsh, versatility is essential, the Tactical Compression Wrap from North American Rescue is the wound-care product preferred by pre-hospital and tactical healthcare professionals. It's an *"any time, anywhere, any application"* *wound-care solution* for medics who want the best there is so they can perform the best job possible.

Features:

- Most versatile pressure bandage system available
- Lightweight and compact, easy to apply
- For extremity, head, shoulder and abdominal wounds, buddy splinting & more
- Tested with great success in Afghanistan and Iraq
- Developed by a Special Forces Medic

Specifications:

- Packaged: L 5.4 in. x W 2.7 in. x D 2 in.
- Bandage: W 5 in. x L 38 in.
- Weight: 6 oz

Item#

30-0016

NSN#

6510-01-614-0938



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS WITH A
MISSION®