

EAGLE IFAK *Individual First Aid Kit*

REV120622

One of NAR's smallest IFAKs, complete with essential life-saving equipment for point-of-wounding and self-aid/buddy-aid care to address penetrating, blast or other traumatic injuries...

A compact, versatile individual first aid kit that offers both MOLLE/PALS – style connectors to attach to the assault vest/body armor and an upper or lower mounted 2.25 in. belt loop attachment to secure the kit in the high or low position on your duty/pistol belt, all in one carrying platform.

Sized as one of the smallest IFAKs in our product line, this innovative kit contains all of the essential equipment for point-of-wounding and self-aid/buddy-aid care to address penetrating, blast or other traumatic injuries. Featuring a Quick Pull Tab, the carrying bag easily opens to reveal the removable kit with just one hand. The kit has a unique, detachable shock-cord inner panel that can be adjusted to carry different components and provides easy access to your IFAK contents by simply pulling the color-coded tab to remove the panel from the carrying bag. The compact size reduces the space requirement needed for an IFAK on your vest or belt.

Features:

- Contains the essential equipment needed for self-aid/buddy-aid to address top leading causes of potentially preventable death in a tactical environment
- Can connect to the assault vest/body armor or to a duty/pistol belt
- Quick Pull Tab easily opens the kit with just one hand
- Has a unique, detachable shock-cord inner panel that can be adjusted to your internal components
- Inner panel is easily removable for easy access to your kit

Specifications:

- Packaged: H 8 in. x W 4 in. x D 3 in.
- Weight: 1 lb 2 oz
- Advanced Weight: 1 lb 3 oz



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

Eagle Individual First Aid Kit

Color	Option	Item #
Black	Basic	85-0177
Coyote	Basic	85-0178
OD Green	Basic	85-0180
MultiCam	Basic	80-1490
Ranger	Basic	80-1491
Black	w/ Combat Gauze®	85-0915
Coyote	w/ Combat Gauze®	85-0982
OD Green	w/ Combat Gauze®	85-1044
MultiCam	w/ Combat Gauze®	85-0191
Ranger	w/ Combat Gauze®	80-1492
Black	Advanced	85-0189
Coyote	Advanced	85-0190
OD Green	Advanced	85-0192
MultiCam	Advanced	80-1495
Ranger	Advanced	80-1496

Available Colors:



PRODUCTS
WITH A
MISSION®

EAGLE IFAK

REV120622

All kits include:



Plus:



Kit Contents:

- 1 x C-A-T® (Combat Application Tourniquet®) Black
- 1 x 4 in. Flat ETD™ Emergency Trauma Dressing
- 1 x HyFin® Vent Chest Seal Twin Pack
- 1 x pair Bear Claw® Nitrile Trauma Gloves, lg.

Basic also includes:

- 1 x Wound Packing Gauze

Combat Gauze also includes:

- 1 x Combat Gauze Hemostatic Dressing

Advanced also includes:

- 1 x Combat Gauze Hemostatic Dressing
- 1 x Nasopharyngeal Airway 28F with Lubricant
- 1 x ARS® Needle Decompression Kit (14 G x 3.25 in.)

The Eagle IFAK's compact size reduces the space requirement needed for an IFAK on your vest or belt.

Quick Pull Tab easily opens the kit with just one hand...



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®