


# NAR CRICOTHYROTOMY KIT *with Bougie Introducer*

REV120822



## ***Pre-inserted bougie for introduction of the tube into the trachea and indication of correct placement prior to insertion...***

The North American Rescue Cricothyrotomy Kit is a convenient device that gives the user a pre-inserted bougie for introduction of the tube into the trachea and indication of correct placement prior to insertion. The bougie-introducer is adjustable for maximum user preference, and the excess can be tucked into the included cap that rests on top of the tracheostomy tube. Contains the same preferred flexible tube as NAR's Tactical CricKit®. This latex free tube is sterile and has a 6.0mm inner diameter and an inflatable cuff.

### **Features:**

- Convenient device that gives the user a pre-inserted bougie for introduction of the tube into the trachea and indication of correct placement prior to insertion
- Adjustable for maximum user preference
- Latex free tube features 6.0mm inner diameter and inflatable cuff
- Contains the same preferred flexible tube as NAR's Tactical CricKit®

### **Specifications:**

- Tube diameter: Inner: 6.0mm | Outer: 8.0mm
- Tube length: 67.0mm
- Packaged:  
L 5.5 in. x H 3.625 in. x D 0.5 in.
- Weight: 0.84 oz

### ***NAR Tracheostomy Kit with Bougie Introducer***

Item#	NSN#
10-0048	



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**