

TACTICAL CRICKIT®

REV080822

The compact all-in-one cricothyroidotomy kit that features all CoTCCC recommendations for a surgical airway kit in the tactical setting



Perform surgical cricothyroidotomy In the tactical environment

Supplying the operator with all the necessary tools to conduct a proper cricothyroidotomy in a low cube and ruggedized package, the Tactical CricKit® is the ideal kit for performing a cricothyroidotomy in austere environments.

For ease in both training and live fire situations, its components are provided in the proper order of use. Once unrolled, the Tactical CricKit's 500-denier pouch provides a clean working surface for the operator, making it viable for use in low-light and battlefield environments.

The Tactical CricKit® also provides the user with a tracheal tube that is both the ideal size and length, along with a securing device to aid in secondary treatment and TACEVAC. To combat microbial contamination, this kit is supplied with alcohol pads for immediate disinfection, and povidone for long-term protection. The Tactical CricKit® was specifically engineered to meet CoTCCC requirements for a surgical airway kit.

Features:

- Provides all of the components necessary to conduct a surgical cricothyroidotomy correctly
- The ideal length and size tracheal tube with securing device included
- Components configured in the sequence of use and ideal for low-light, battlefield environments

Specifications:

- Kit Unrolled -
H 6.5 in. x W 7.5 in. x D 1.25 in.
- Packaged -
H 6.5 in. x W 2.5 in. x D 1.25 in.
- Weight: 3.7 oz

COY

**Available
Color:**

Tactical CricKit®



Color	Item#	NSN#
COY	10-0017	6515-01-540-7568



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

TACTICAL CRICKIT®

REV080822



Kit Contents:

- 1 x Rolled Nylon Pouch
- 1 x Alcohol packet
- 1 x Povidone Iodine packet
- 1 x Gauze Pad (4 in. x 4 in.)
- 1 x Protected #10 Scalpel
- 1 x NAR Tracheal Hook
- 1 x 10cc Syringe
- 1 x Cric. Tube with Obturator (Flexible, Cuffed, 6mm)
- 1 x Tube Securing Strap
- 1 x Plastic Vacuum Bag



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®