

CELOX™ FIRST AID Z-FOLD HEMOSTATIC DRESSING

REV011223

Solutions For First Aid

Stop bleeding from minor cuts, lacerations and abrasions at home, at work, for leisure, for travel. Be prepared with CELOX™ First Aid Hemostatic Gauze Products. CELOX™ Granules are bonded to the surface of CELOX™ Gauze, which has a more stable, denser consistency than traditional thready gauze and does not compact under pressure. Just pack or cover the bleeding area and apply pressure.

At Home Or On The Go

Accidents always happen at the worst times but Celox™ First Aid can help to get control of bleeding injuries. Lightweight, compact and affordable – perfect for the kitchen, at the workbench, sports, hobbies, camping and whenever you are on the move.

Easy To Use

CELOX™ First Aid Hemostatic Gauze is now Z-Folded for easy use – tear it open, apply to the source of bleeding, cover and apply pressure. It's that simple!

Features:

- Advanced hemostatic agent combined with the simplicity of gauze
- Fast-acting technology
- Easy to apply on superficial wounds, minor cuts, and abrasions
- 3 in. x 5 ft. Z-Fold configuration

Dimensions:

- Packaged - L 5.5 in. x W 4.5 in.
- Gauze - L 5 ft. x W 3 in. (Z-Fold Configuration)
- Weight: 0.10 lbs

CELOX™ First Aid Z-Fold Hemostatic Dressing

Item#	NSN#
30-0300	---



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS

WITH

MISSION®