

HYFIN® CHEST SEAL *Occlusive Dressing*

HyFin® Chest Seal creates a superior seal...

The new HyFin® Chest Seal is a completely occlusive dressing used for the treatment of penetrating injuries to the chest. Now featuring a larger, 6 in. x 6 in. size and enhanced adhesive gel formulation, the new HyFin® Chest Seal creates a superior seal, even in the most adverse conditions. The new HyFin® Chest Seal also features a transparent backing to allow placement directly over the wound and easily conforms to the casualty's chest. It also has a large, easily-visible red elliptical tab for single step peel-and-apply application that can also facilitate the burping of a wound if necessary.

The HyFin® is Made in America and is 100% latex free. To ensure high quality performance, NAR packages each chest seal using its' signature Red Tip Technology® red tear notches in a rugged, easy-to-open, foil pouch. Maximum success in managing open chest wounds requires an effective, easy-to-use and reliable chest seal. The HyFin® Chest Seal addresses each of these requirements and more.

Features:

- New, more compact packaging
- Advanced adhesive technology
- Oversized Red elliptical tab for single step peel-and-apply application that can also facilitate burping of the wound
- Conforms to the casualty's chest
- Effective in extreme environments
- Seal can be written on for patient treatment documentation if necessary
- Latex-free
- Made in America

Dimensions:

- Chest Seal: H - 6.0 in. x W 6.0 in. x D 0.1 in.
- Packaged: H 4.625 in. x W 7.5 in. x D 0.1 in.
- Weight: 1.44 oz



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277



HyFin Chest Seal Individual Occlusive

Item #	NSN#
10-0015	6510-01-532-8019

Unique adhesive gel formulation does not damage mannequins or leave behind excessive adhesive residue

PRODUCTS WITH A
MISSION®