

# CRISIS INCIDENT RESPONSE KITS

REV120618



ORG BLK



***Designed to treat penetrating trauma and major bleeding as well as rescue breathing and basic first aid***

The Crisis Incident Response Kit from North American Rescue is designed to provide the necessary medical products to treat the more serious medical emergencies that may occur in a school, public or large setting.

This ruggedized Nylon kit was designed to combat life threatening emergencies such as the inability to breathe or major bleeding.

This rugged nylon kit provides carrying options such as over the shoulder, hand carry and fanny-pack style.

### Features:

- The necessary emergency first aid products to treat several injured persons ranging from minor to severe injuries
- The School First Aid Kit designed to treat penetrating trauma and major bleeding as well as the administration of rescue breathing & basic first aid
- Rugged carrying bag constructed of durable nylon
- Multiple carry options: over the shoulder, hand carry & fanny pack style

### Dimensions:

- Packaged:  
H 7.75 in.  
W 13 in.  
D 9 in.
- Weight: 5 lb

### Crisis Incident Response Kit

Color	Item #	NSN#
ORG	85-0410	
BLK	85-1338	



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**

# CRISIS INCIDENT RESPONSE KITS

REV120618



## Kit Contents:

- 2 x C-A-T® Tourniquet -  
ORG TQ in #85-0410 | BLK TQ in #85-1338
- 4 x NAR Compressed Gauze
- 4 x NAR 4 in. Flat (ETD)  
Emergency Trauma Dressing
- 4 x NAR 6 in. (ETD)  
Emergency Trauma Dressing
- 1 x NAR Abdominal (ETD)  
Emergency Trauma Dressing
- 2 x Hyfin® Vent Compact Twin Pack
- 2 x CPR Microshield
- 2 x SAM Splint II
- 1 x roll 2 in. Surgical Tape
- 5 x pair Bear Claw Gloves, large
- 1 x NAR 7.25 in. Trauma Shears
- 50 x Flexible Fabric Bandages
- 2 x Survival Blanket

*Ruggedized Nylon kit  
designed to address life  
threatening emergencies  
such as major bleeding or  
the inability to breathe*



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**