

CASUALTY RESTRAINT STRAPS

TACEVAC

2 pc Plastic Buckle

Color	Item#	NSN#
BLK	50-0023	
ORG	50-0024	

Dimensions:

L 5 ft. x W 2 in.

Wt: 4 oz



2 pc Metal Buckle/Speed Clip

Color	Item#	NSN#
BLK	50-0029	
ORG	50-0030	

Dimensions:

L 5 ft. x W 2 in.

Wt: 11.7 oz



2 pc Plastic Buckle/Speed Clip

Color	Item#	NSN#
BLK	50-0026	
ORG	50-0025	

Dimensions:

L 5 ft. x W 2 in.

Wt: 7.8 oz



2 pc Metal Buckle

Color	Item#	NSN#
BLK	50-0027	
ORG	50-0028	

Dimensions:

L 5 ft. x W 2 in.

Wt: 8 oz



North American Rescue's Casualty Restraint Straps are constructed using heavy-duty nylon webbing and convenient seat belt style buckles which provide the utility you need for quick and secure immobilization with any standard back board or basket stretcher.

PRODUCTS WITH A
MISSION[®]



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

CE

BLK ORG