

# ANKLE TOURNIQUET HOLSTER



## Lightweight and Comfortable...

The Ankle Tourniquet Holster (ATH) is designed to provide the capability for every law enforcement officer or other individuals who operate in potential high-threat environments to easily carry their C-A-T® Tourniquet on their person, where it is needed most.

Made of 3mm perforated neoprene and weighing only 2.5 oz, the ATH is a lightweight and comfortable tourniquet carrying device that allows easy access and quick deployment of the life-saving C-A-T® Tourniquet.

It can be worn directly over your high top duty boot or on your ankle. The ATH is easily placed via hook and loop attachment and stays in-place during the performance of your duties. The ATH firmly holds your C-A-T® in the holster while allowing rapid access with either hand for fast deployment when needed.

## Features:

- Allows easy access and quick deployment of tourniquet
- Made of 3mm perforated neoprene and weighing 2.5 ounces
- Hook and loop attachment ensure that Holster stays in-place
- Allows rapid access with either hand for fast deployment when needed
- Accommodates both Gen 6 and Gen 7 C-A-T® Tourniquets (not included)

## Dimensions:

- Size: L 16 in. x W 4 in. x D 0.25 in.  
(One Size Fits Most)
- Weight: 2.5 oz

NOTE: Tourniquets sold separately



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277