

TRUECLOT® Hemorrhage Training Kit™

REV072318

The TrueClot® Hemorrhage Training Kit™ by Luna is designed with everything you need to demonstrate and train on advanced wound packing techniques with the use of advanced hemostatic gauzes. This highly realistic Wound Packing Trainer features a large caliber gunshot wound with hidden bleed points and a simulated fracture/defect hard plastic bone insert to aid in teaching deep packing techniques necessary for larger wounds, ie. “pack to the bone.” As Luna’s proprietary TrueClot® Clotting Blood Simulant is pumped into the wound track using the attached 1 liter squeeze bottle and makes contact with the Simulated Hemostatic Gauze, realistic simulated blood clots form on the gauze at the point of bleeding to demonstrate the reaction of Hemostatic Gauze in a wound. Mimicing the mechanics of actual wound packing, the trainer can be fully packed with 12 ft of gauze and have external pressure applied to help slow blood flow.

Features:

- Durable, highly realistic Wound Packing Trainer Kit
- Task Trainer features a large caliber gunshot wound with hidden bleed points and simulated fracture/defect hard plastic bone insert to aid in teaching deep wound packing techniques
- Includes TrueClot® Clotting Blood Simulant and Simulated Hemostatic Gauze which demonstrates the reaction of Hemostatic Gauze in a wound
- Durable design features a sturdy silicone construction with no moving parts, reinforced tubing mount, and a more realistic wound pattern and outer anatomy – ideal for the repeated use of a training environment

Dimensions:

- H 7 in. x W 5 in. x D 2.5 in.
- Weight: 7 oz

TrueClot® Hemorrhage Training Kit™

Item #	NSN#
83-0020	

Introduce students to realistic hemorrhage control for **TCCC, TECC, TEMS, PHTLS, Active Shooter Training** and **general wound packing techniques** for **advanced hemostatic gauze training**



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**

TRUECLOT® Hemorrhage Training Kit™

REV072318

 <p>QTY 1 Lg. GSW Task Trainer</p>	 <p>QTY 4 TrueClot® Blood Simulant</p>	 <p>QTY 9 TrueClot® Hemostatic, 12 ft</p>
 <p>QTY 1 Bottle, Delivery with Tube & 1-way Valve</p>	 <p>QTY 1 Bottle, Mixing & Storage</p>	 <p>QTY 1 Hard Case with Instructions</p>

Ideal for training wound packing techniques featuring a durable silicone construction with no moving parts, reinforced tubing mount, and a more realistic wound pattern and outer anatomy...



Kit Contents:

- 1 x Large Caliber Gunshot Wound Task Trainer with Bone Insert
- 9 x TrueClot® Simulated Hemostatic Gauze Pack, 12 ft
- 4 x TrueClot® Blood Simulant Concentrate Pack, 1 liter (when mixed)
- 1 x Blood Simulant Delivery Bottle & Tubing with one way valve
- 1 x Blood Simulant Mixing/Storage Bottle
- 1 x Instructions for Use
- 1 x Rigid Carrying Case



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A 
MISSION®