

# FIBER OPTIC LARYNGOSCOPE BLADES



◀ Mac #2



◀ Mac #3



◀ Mac #4

## Fiber Optic Laryngoscope Blades, Mac:

For improved safety and convenience, use the Mac Fiber Optic Laryngoscope Blades from North American Rescue. These blades offer the superior illumination of fiber optic technology and compatibility with any standard (green spec) fiber optic handle. But, that's not all. Enhanced polymer construction gives them exceptional strength and rigidity. More importantly, this construction reduces the risk of damage to the patient's teeth or soft tissue.

### Features:

- Single-use laryngoscope blades compatible with any standard fiber optic handle
- Enhanced polymer construction for less risk of damage to patient's teeth or soft tissues
- Exceptional strength and rigidity plus superior illumination
- Designed for safety and convenience

### Mac #2:

- L 4.6 in. x W 1.4 in. x D 1.1 in.
- Weight: 0.8 oz

### Mac #3:

- L 5.25 in. x W 1.5 in. x D 1 in.
- Weight: 1 oz

### Mac #4:

- L 6 in. x W 2.5 in. x D 1 in.
- Weight: 1 oz



Size	Item#	SEL#	NSN#
Mac #2	10-0012	09ME-02-AWMG	6515-01-460-4686
Mac #3	10-0013	09ME-02-AWMG	6515-01-460-4685
Mac #4	10-0014	09ME-02-AWMG	6515-01-460-4672



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

PRODUCTS WITH A  
**MISSION®**