

HELO TAGLINE KIT

REV071520

Help control litter movement during hoist operations – a requirement on all Army MEDEVAC helicopters

The Helo Tagline Kit from North American Rescue contains 250 ft of 7.5 mm Triton Pegasus® Rope featuring a highly visible polypropylene sheath with a super strong Innegra™ core.

This kit includes 1 locking carabiner, 10 weak links (designed to break away if the litter is snagged during the hoist) and 1 stainless steel quick link..

Features:

- Designed to help control litter movement during hoist operations
- Contains 250 ft of 7.5 mm Triton Pegasus® Rope
- Includes 1 locking carabiner, 10 weak links (designed to break away if the litter is snagged during the hoist) and 1 stainless steel quick link

Kit Contents:

- 250 ft of 7.5 mm Triton Pegasus® Rope
- 1 x rockD Lock Carabiner
- 1 x Stainless Steel Quick Link
- 10 x Weak Link

Dimensions:

- L 11 in. x W 7 in. x D 7 in.
- Weight: 5 lb 13 oz

Helio Tagline Kit

Item#	NSN#
90-0069	

Helio Tagline Kit Contents:



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH  **MISSION®**