


# T2 TACTICAL TRIAGE KIT

Low-cube, vacuum-sealed kit to document the initial triage of up to 4 casualties



## Ultimate Source for On-Scene Triage...

A compact incident command system, the Tactical Triage Kit streamlines medical triage. A triage flowchart and casualty counter are laminated back-to-back and attached to a break-away lanyard for peak efficiency.

The flowchart summarizes the SALT triage process in easy-to-follow steps. The background of the casualty-count card is divided into colored sections matching the casualty-status levels defined by SALT (IMMEDIATE—red, DELAYED—yellow, MINIMAL—green, EXPECTANT—gray, and DEAD—black). A black grease pencil lets users quickly check off numbers in each category.

All this critical triage equipment is contained in a lightweight, rugged vacuum-sealed pouch, making it ideal for trauma packs and vehicle kits.

### Features:

- SEL#: 09ME-01-MC1K
- Compact vacuum-sealed pouch (10.5 in. x 7 in.)
- Laminated triage flow sheet & casualty count card w/break-away lanyard
- 4 NAR T2 Individual Kits with liquid-resistant triage/treatment tag, plastic sheath, chemical light, pencil, and status stickers

### Kit Contents:

- 1 x Triage Flow Sheet and Casualty Count Card (laminated) with break-away lanyard
- 4 x NAR T2 Tag®
- 1 x Grease Pencil (Black)

### Dimensions:

- Packaged - L 5.25 in. x W 7 in. x D 1.25 in.
- Weight: 4.9 oz



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277