

# NAR RESCUE SLEEVE II™ & Accessories

REV090120



NOTE: Spine board sold separately.  
Hoisting straps & carabiners not included.

BLK  
ORG

SEL#: 09ME-04-SPIN  
SEL#: 03WA-02-BORD

## The Ultimate In Patient Immobilization...

The NAR Rescue Sleeve II™ is the industry's ultimate immobilization system supporting both vertical and horizontal hoisting. Simple and secure, the Rescue Sleeve II™ gives confidence to operators and comfort to patients. Used extensively in the field and refined for enhanced stability and overall functionality, it is standard issue to all U.S. Navy ships.

### Features:

- Vertical and horizontal lifting-strap attachments
- Removable head-securing blocks
- Adjustable straps secure head and chin
- Hoist-cable storage space
- Padded head flaps
- Backboard compartment
- Padded carrying handles
- Easy to apply
- Replaceable parts are also available:  
**Head Block** (#50-0042)  
**Head and Chin Straps** (#50-0001 | NSN# 6515-01-477-8921)  
**Carrying Case** (# 50-0002 | NSN# 4240-01-477-9408)
- Standard issue to all U.S. Navy ships



### Specifications:

- Folded:  
L 25.5 in. x W 13 in. x D 8 in.
- Open:  
L 73 in. x W 18 in. x D 2 in.
- Payload: 1000 lb (453.6 kg)
- Weight: 16 lb

### NAR Rescue Sleeve II

| COLOR  | ITEM#   | NSN#             |
|--------|---------|------------------|
| Black  | 50-0011 | 6530-01-531-6649 |
| Orange | 50-0012 | 6530-01-477-8525 |



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**