

# MEDIC LEG RIG™ Kit

REV092121

## Designed for the Tactical Medic

Tactical Healthcare Professionals who want rapid access to critical basic and advanced equipment by their side choose the Combat Casualty Response Kit®—Medic/Leg Rig™. The container/bag is constructed of 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction, Ghillie-Tex IR signature reduction hardware and silent, self repairing nylon zippers with 550 pull cords.

The kit can be employed in three ways

1. Using the split-tale connector with height adjustment and Quick-Grip™ leg straps
2. PALS/MOLLE style webbing for att. to the aidbag/ backpack
3. Fasten to the web/gun belt using the quick weave method.

Belt and leg attachment straps can be removed and reattached as needed to allow maximum flexibility during load carriage configuration.

### Features:

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- Ghillie-Tex IR signature reduction hardware
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE Style Attachment System Compatible
- Trauma Shear Storage on Back for Quick Access

### Specifications:

- Case Closed: H 9.5 in. x W 9 in. x D 7.75 in.
- Weight: 3 lb 9 oz



### Available Colors:



Color	Item#	NSN#
BLK	80-0046	6545-01-543-3185
ODG	80-0047	6545-01-543-3174
COY	80-0048	6545-01-543-3170
DUC	80-0087	
MTC	80-0249	



**NORTH AMERICAN RESCUE®**


[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**

# MEDIC LEG RIG™ Kit

REV092121

*Rapid access & organization of critical equipment for basic to advanced medical field care*



## Kit Contents:

- 5 x pair Black Talon® Nitrile Gloves
- 2 x Nasopharyngeal Airway 28F with Lubricant
- 1 x Petrolatum Gauze (3 in. x 18 in.)
- 1 x HyFin® Chest Seal, Occlusive Single
- 2 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 2 x C-A-T® (Combat Application Tourniquet®)
- 2 x ETD™ 6 in. Emergency Trauma Dressing
- 4 x S-Rolled Gauze™ (4.5 in. x 4.1 yd)
- 2 x NAR Saline Lock Kit
- 1 x BOA® IV Constricting Band
- 1 x Trauma Shears (7.25 in.)
- 4 x Polycarbonate Eye Shield
- 1 x Scissor Leash™
- 1 x Surgical Tape (2 in.)
- 1 x Combat Casualty Reference Card
- 2 x Triage Card (NAR T2 Tag®)



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**