


MASSIVE HEMORRHAGE | AIRWAY | RESPIRATION | CIRCULATION | HYPOTHERMIA

# EMS *Catalog*





**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

Since 1996, North American Rescue has focused on **increasing survivability** on both the battlefield and the homefront.

With our innovative medical/rescue equipment, NAR continues to lead the way by delivering the **tailored, mission-critical healthcare solutions** needed by today's first responders to counter the rapidly evolving threats of the 21st century. See what a lifesaving difference our combat proven **Products with a Mission®** make in the prehospital environment!





We are dedicated to decreasing preventable death by providing our products at special pricing to our EMS and pre-hospital life savers. To see if you qualify for this discount, please contact us or fill out our online NAR Pricing Application. Our customer service department is standing by to help!

# M.A.R.C.H.

Massive Hemorrhage ..... **4-6**

Airway Management ..... **7**

Respiratory Support ..... **8-10**

Circulation & Vascular Access ..... **10-13**

Hypothermia Management ..... **14**

## Trauma Kits

IFAKS ..... **14-18**

Custom Kitting Solutions ..... **15**

Multiple Injury Kits ..... **14,18**

## Extras

Splinting & Evacuation ..... **19**

Treatment Support Accessories ..... **21-22**

Education & Training ..... **20,23**

# Massive Hemorrhage | Tourniquets

## C-A-T® Combat Application Tourniquet Gen 7

- **Single route buckle** decreases application time to help decrease blood loss, allows for slack removal for proper TQ application
- **Applied same way** to both upper and lower extremities
- True **one-handed tourniquet**
- **Scientifically proven** 100% effective in occluding blood flow in both upper and lower extremities by the USAISR and recommended by CoTCCC



**The FASTEST, SAFEST,  
Most EFFECTIVE  
Prehospital Field Tourniquet!**

### SPECIFICATIONS

Packaged: H 6.5 in. x W 2.4 in. x D 1.5 in.  
Device: L 37.5 in.  
Weight: 2.7 oz



**online video available**

< scan to view

### COMBAT APPLICATION TOURNIQUET

COLOR	STYLE	ITEM #	NSN #
● Black	Tactical	30-0001	6515-01-521-7976
● Orange	Rescue	30-0023	6515-01-612-7129
● Blue	Trainer	30-0033	6910-01-560-2972

## RIGID C-A-T® TQ CASE with Shirt Shield

- Allows you to carry a Gen 7 C-A-T® Tourniquet on your duty, patrol, or pistol belt so you have life-saving equipment on your person at all times
- Includes a Shirt Shield complete with mounting hardware to protect your uniform shirt from undue wear caused by the friction from the windlass clip
- Blade-Tech TekLok Belt Attachment fits up to a 2.25 in. belt
- Can be mounted either vertically or horizontally or attached to a Safariland Drop Leg Shroud
- Injection Molded Nylon Case stands up to rugged LE/MIL usage
- Cross medical identifier patch included



NOTE: All Tourniquets sold separately

### RIGID C-A-T TQ CASE with Shirt Shield

COLOR	ITEM #
● Black	80-0944
● Black, Basket Weave	80-0962
● Black, Polished	80-0963
● Coyote	80-0945
● Ranger Green	80-0946
● Navy	80-1161

\* More colors available at [NARescue.com](http://NARescue.com)

### SPECIFICATIONS

H 5.88 in. x W 2.5 in. x D 2 in. | Weight: 5.36 oz

# LEATHER C-A-T® HOLDERS

- Designed to meet the specific duty gear needs of law enforcement officers - for use on daily duty belt to fit within organizational uniform standards
- Securely fastens to your duty belt with easy one-handed access to allow you to take control of extreme medical situations in seconds
- Constructed of black full grain leather in different finishes, attachment options, and snaps colors
- High quality snaps ensure secure placement on your duty belt and easy accessibility to your tourniquet
- Ideal for the C-A-T® Gen 7 but also accommodates other popular windlass TQs
- Attachment options: Clip, Loop or Hook & Loop Strap attachment

## SPECIFICATIONS

L 5.75 in. x W 2 in. x D 2 in.  
Weight: 3.2 oz

## UNIFORM DUTY BELT LEATHER TQ HOLDER

FINISH	ATTACHMENT	SNAP	ITEM #
Basket Weave	Clip	Silver	30-0192
High Gloss	Clip	Silver	30-0194
Full Grain	Clip	Silver	30-0193
Basket Weave	Loop	Silver	30-0195
High Gloss	Loop	Silver	30-0198
Full Grain	Loop	Silver	30-0197
Basket Weave	Strap Attach	Silver	30-0228
High Gloss	Strap Attach	Silver	30-0229
Full Grain	Strap Attach	Silver	30-0196



See website for more  
Uniform Duty Belt Tourniquet Holder options

## JETT® Junctional Emergency Treatment Tool

- Immediately treat high extremity wounds in the thigh/groin area that is not amenable to TQ application
- Simultaneous occlusion of blood flow to both lower limbs that will not impede respiration
- Patent 9,492,177

## SPECIFICATIONS

Packaged:  
H 5.5 in. x W 7 in. x D 3.75 in.  
Overall Device: L 56.75 in.  
Buttocks Pad: L 19 in. x W 6.5 in.  
Compression Pad:  
L 4.85 in. x W 2.5 in. x D 11 in.  
T-Handle:  
H 4.62 in. x W 4.04 in. x D 0.63 in.  
Weight: 1 lb 9.6 oz



## JETT

ITEM #	NSN #
30-0088	6515-01-616-5841

online video  
available  
scan to view >



# Massive Hemorrhage | Bandages & Gauze

## NAR COMPRESSED GAUZE

- Low-cube, ultra-compact sterile gauze
- Premium rolled 100% cotton gauze for maximum absorption



### COMPRESSED GAUZE

STYLE	ITEM #	NSN #
Responder	30-0182	----
Standard	30-0052	6510-01-503-2117

### SPECS

Gauze: W 4.5 in. x L 4.1 yd  
 Packaged: L 3 in. x W 2.5 in. x D 1 in.  
 Weight: 1.2 oz

## NAR WOUND PACKING GAUZE

- Ultra compressed, Z-fold gauze indicated for wound packing
- Soft, white 1-ply gauze with an embedded x-ray detectable line



### SPECS

Gauze: W 3 in. x L 5 yd  
 Packaged:  
 H 3 in. x W 2 in. x D 1 in.  
 Weight: 12 oz

### WOUND PACKING GAUZE

ITEM #	NSN #
30-0054	6510-01-695-5824

## RESPONDER ETD™ Emergency Trauma Dressing

- Sterile pressure bandages that deliver extremely effective pressure & excellent absorption
- Suitable for applying immediate direct pressure
- Equipped with a non-adherent pad
- Easy-to-use hook & loop securing device with no confusing pressure bars or hooks
- Simple application with the Quick grip roll control system which prevents the dressing from unrolling if dropped
- Rugged vacuum-sealed, low-cube packaging with easy-open tear notches



### RESPONDER ETDs

STYLE	BANDAGE	PAD	WGT	PACKAGE	ITEM #
4 in. Flat	W 4 in. x L 60 in.	4x6	2 oz	L 4.25 in. x W 2.75 in. x D 0.88 in.	30-0171
4 in. Flat (Mini)	W 4 in. x L 42 in.	4x6	2.4 oz	L 5 in. x W 3.5 in. x D 0.75 in.	30-0179
6 in. Rolled	W 6 in. x L 60 in.	6x6	2.9 oz	L 6 in. x W 2.5 in. x D 2.1 in.	30-0168

# Massive Hemorrhage | Hemostatics

## COMBAT GAUZE® LE | COMBAT GAUZE®

- Sterile, non-woven Z-fold gauze impregnated with kaolin for temporary external control of bleeding
- Effective, simple and safe to use on wounds
- Features an embedded x-ray detectable strip
- Does not generate heat



### COMBAT GAUZE Z-FOLD

VERSION	ITEM #	NSN #
Combat Gauze	30-0039	6010-01-562-3325
Combat Gauze XL	30-0070	6510-01-600-4705
Combat Gauze LE	30-0109	----

### SPECIFICATIONS

#### Standard & LE Gauze:

L 4 yd x W 3 in.  
 Pkg: H 5 in. x W 5.5 in. x D 0.625 in.  
 Weight: 0.8 oz

#### XL Gauze:

L 4 yd x W 4 in. (2 ply)  
 Pkg: H 6.5 in. x W 6.5 in. x D 0.75 in.  
 Weight: 1.75 oz

## BAC-PACK™ *Bougie-Aided Cricothyroidotomy Pack*



- All-in-one cricothyroidotomy kit for use in any environment
- Pre-inserted bougie in tube allows for a more rapid, easy procedure
- Tube meets TCCC recommendations for bougie-aided open surgical technique
- #10 Scalpel included for safety and ease of use
- Small, durable pouch keeps the items organized and ready to use
- Lightweight and ruggedized
- Components configured by order of use
- Scalpel, tube, syringe all provided sterile
- Patent pending

### SPECIFICATIONS


Packaged: L 5.5 in. x H 3.625 in. x D 0.5 in.  
 Inner tube diameter: 6.0mm  
 Outer tube diameter: 8.0mm  
 Tube length: 67.0mm  
 Weight: 0.84 oz  
 US Patents: Pending

### BAC-PACK



ITEM #

10-0052

### INCLUDED IN THE BAC-PACK:



## NAR Tracheostomy Kit with BOUGIE INTRODUCER



### TRACHEAL TUBE with FLANGE

### SPECIFICATIONS


Tube diameter:  
 Inner: 6.0mm | Outer: 8.0mm  
 Tube length: 67.0mm  
 Packaged:  
 L 5.5 in. x H 3.625 in. x D 0.5 in.  
 Weight: 0.84 oz

- Convenient device that gives the user a pre-inserted bougie for introduction of the tube into the trachea and indication of correct placement prior to insertion
- Adjustable for maximum user preference
- Latex free tube features 6.0mm inner diameter and inflatable cuff
- Contains the same preferred flexible tube as NAR's Tactical CricKit®

### TRACHEOSTOMY KIT

ITEM #

10-0048




## HYFIN® CHEST SEALS



### HYFIN® VENT COMPACT, *Twin Pack*

- Two Chest Seals packaged approximately 25% smaller than our standard HyFin®
- 3-channel pressure relief vents
- Two Chest Seals designed for the treatment of entry/exit or multiple penetrating wounds
- Large red tab designed for single step, peel and apply application and for burping the wound if necessary



#### HYFIN VENT COMPACT

ITEM #
10-0042

#### SPECIFICATIONS

Device (excludes tab): 4.75 in. x 4.75 in.  
 Weight: 155 oz  
 Packaged: H 3.88 in. x W 6.6 in. x D 0.13 in.


online video available

scan to view >



### HYFIN® VENT *Twin Pack*

- Two individually packaged, vented/valve-type chest seals in one unique package designed for treatment of both entry/exit wounds or multiple penetrating injuries
- Patented 3 channel pressure relief vents
- Ruggedized dual perforated pouch designed to allow the provider to open only one dressing at a time, as needed
- Advanced adhesive technology designed to provide a superior seal



#### SPECIFICATIONS

Device (excludes tab): 6 in. x 6 in. Wgt: 2.5 oz  
 Packaged: H 7.5 in. x W 4.5 in. x D 0.25 in.  
 Unfolded: H 7.5 in. x W 9 in. x D 0.13 in.

#### HYFIN VENT TWIN PACK

ITEM #	NSN #
10-0037	6510-01-642-6212



Includes 6 ft 6 in. of oxygen tubing & rugged, hard plastic, twist-top case

### CYCLONE® Improved **POCKET BVM**

- Unique collapsible design packaged 25% smaller than conventional BVMs
- Quick to deploy and does not need to be manually inflated prior to use

#### SPECIFICATIONS

Diameter: 5.3 in. x H 2.9 in.  
 Weight: 17.4 oz


#### CYCLONE® **POCKET BVM**

ITEM #	NSN #
10-0019	6515-01-568-0193



## 10 ga x 3.75 inch **SPEAR™** Simplified Pneumothorax Emergency Air Release

- Latest innovation in Needle Decompression technology
- Designed based on known gaps in the management of tension pneumothorax
- Enhanced needle tip for improved insertion control
- User friendly features like spin-lock connection, depth markings, and detachable low-pressure one-way valve
- Fenestrations along catheter for additional air flow
- Familiar, rugged, pen-style case for ease of storage and deployment
- Allows for either lateral decompression or anterior needle thoracostomy
- Catheter made of flexible kink-resistant material with a soft tip
- Single use only. Supplied sterile



### SPEAR

GAUGE	ITEM #	NSN #
10 ga	10-0051	6515-01-699-4360



online video available  
scan to view >



## 3.25 inch **ARS®** NEEDLE DECOMPRESSION KIT

- Penetrate the pleural space in 99% of casualties presenting the signs and symptoms of a tension pneumothorax
- Quality, reliable devices with the appropriate length needle/catheter
- Capless flash chamber provides immediate confirmation of correct placement
- Patented design packaged in an easy-to-open, crush-proof protective tube featuring easy-to-identify textured twist top with clip



### ARS NEEDLE DECOMPRESSION KIT

GAUGE	ITEM #	NSN #
14 ga	ZZ-0056	6515-01-541-0635
10 ga	ZZ-0298	6515-01-673-1701

### SPECIFICATIONS

Pkg: H 5.7 in. x W 0.7 in. x D 0.8 in.  
Weight: 0.6 oz



online video available  
< scan to view

U.S. Design Patents: D584,409 S; D595,847S; D584,410S; U.S. Patent 7,874,426B2; © North American Rescue,® LLC., 127820, 127822, 127823, 001002372-0001, 001013940-0001, 001013940-0002



## NAR SIMPLE THORACOSTOMY KIT


- Compact kit that contains all necessary equipment to perform simple thoracostomy
- Includes 3-pack skin prep for performing the procedure bilaterally if needed
- HyFin Compact Vented Twin Pack Chest Seal to cover site if required

### SIMPLE THORACOSTOMY KIT

ITEM #	UNIT
83-0027	Single

### SPECIFICATIONS

Packaged: L 4.5 in. x W 7.5 in. x D 0.75 in.  
Weight: 5.28 oz



QTY 1 Rochester Pear Curved 6.25" Hemostat



QTY 1 Protected #10 Scalpel



QTY 2 Gauze Pads 4x4" (2pk)



QTY 1 HyFin® Vent Compact, Twin



QTY 1 Chloraprep Swabsticks, 3pk



QTY 1 Permanent Marker, Small

## CHEST TUBE INSERTION KIT

- Designed to provide the necessary equipment needed to perform emergency chest tube thoracostomy
- Compact and lightweight kit packaged in an easy to open, vacuum sealed pouch
- Can be stored in your kit or crash cart



### CHEST TUBE INSERTION KIT

ITEM #	UNIT
83-0007	Single

### SPECIFICATIONS


Packaged: L 10 in. x W 6 in. x D 2.25 in.  
Weight: 11.2 oz



QTY 1 Gauze, Petrolatum



QTY 1 Protected #10 Scalpel



QTY 1 5.5" Kelly Hemostat, Cvd



QTY 1 Drape, Fenestrated



QTY 1 Surgical Gloves, Latex



QTY 1 Valve Chest Tube



QTY 1 36 French Chest Tube




QTY 1 Chloraprep Swabsticks, 3pk



QTY 1 4" x 4" Split Sponge



QTY 1 O-Silk Sutures, 30 in.



QTY 1 8" Curved Roch-Peon

## Circulation & Vascular Access

### BOA® IV Constricting Band

- Direct peripheral IV access utilizing a revolutionary rolling action for transfusing venous blood to the distal veins while still allowing arterial flow
- Quick Click connector & two-finger Quick Release improves "time on task"
- XL version stretches approximately 50% longer than the standard BOA®

### SPECIFICATIONS


Standard: Device - L 7.5 in. to 18 in. Weight: 1.2 oz	XL: Device - L 10 in. to 27 in. Weight: 1.5 oz	Packaged: L 4 in. x W 3 in. x D 0.9 in. US Patents: 5,304,202 & 5,607,448
-------------------------------------------------------------	------------------------------------------------------	---------------------------------------------------------------------------------



### BOA IV CONSTRICTING BAND

SIZE	ITEM #	NSN #
Standard	30-0009	6515-01-537-2611
XL	30-0071	----

## SAM<sup>®</sup> IO Intraosseous Access System



### INSERTION SITES:

- Proximal Humerus (left or right)
- Proximal Tibia (left or right)
- Distal Tibia (left or right)
- Distal Femur (pediatric patients)

## SAM<sup>®</sup> IO KIT *Intraosseous Access*

- Designed for intuitive, effective IO placement
- Includes multi-use manually operated intraosseous access driver with passive ratcheting mechanism that rotates clockwise when resistance is encountered
- Provides controlled insertion of the IO needle with a tactile feel of cortical penetration
- Includes 15 mm, 25 mm, and 45 mm Needle sets containing extension tubes and single-use sharps containment to address a broad range of patient demographics and multiple application sites
- Featuring a stabilizer with cylinder design intended to protect the catheter from dislodgement

### SPECIFICATIONS

H 7 in. x W 11.25 in. x D 5.5 in.  
Weight: 2.32 lb

### SAM IO KIT

#### ITEM #

30-0307


## SAM<sup>®</sup> MEDICAL

See item# below for additional **SAM IO System** components also available for purchase:

### SAM IO INTRAOSSEOUS ACCESS SYSTEM / INDIVIDUAL ITEMS

PRODUCT	UNIT	ITEM #
Driver	Single	30-0255
Needle 15mm	Single	30-0257
Needle 15mm	5-Pack	30-0258
Needle 25mm	Single	30-0259
Needle 25mm	5-Pack	30-0260
Needle 45mm	Single	30-0261
Needle 45mm	5-Pack	30-0262
Stabilizer	Single	30-0263
Stabilizer	5-Pack	30-0264
Adaptor*	Single	30-0253
Training Kit	Kit	30-0256
Storage Case	Single	30-0254


\*SAM IO Adaptor is designed to accommodate alternative IO needle designs with the SAM IO Driver



### SAM IO KIT



Includes **Custom Field Storage Case** with elastic loops for easy transport



Extension Tube & NeedleVISE<sup>®</sup> included in each needle set

## LIFEWARMER™ Quantum™ IV Fluid & Blood Warming System

- **Lightweight** The Quantum™ has the lightest operational weight of any portable blood & fluid warming system (Total Operational Weight: 22 oz)
- **Powerful** Warms two units of 39°F (4°C) refrigerated whole blood at 100mL/min+ on a single charge
- **Compact** Fewest components to carry (three)
- **Fast** Warm-up time of only 24 seconds
- **Intuitive** The thermal transfusion/infusion set looks, handles, and sets up/primes like any conventional IV set
- **Safe** Tested safe for blood, blood products & IV solutions

### SPECIFICATIONS

Total Operational Weight: 22 oz  
(Quantum™ Controller, Battery & TTS-B)

**Battery:** 16 oz  
**Controller:** 3 oz  
**Battery Charger:** 22 oz  
**Controller Cover:** 0.23 oz  
**Thermal Transfusion Set (TTS-B):** 2.7 oz

### COMPLETE QUANTUM IV & BLOOD WARMING SYSTEM CONTENTS:



### BLOOD or BLOOD PRODUCTS

From a **single battery charge**, capable of warming **1000 mL** of **39°F (4°C) Whole Blood** to an output temperature of **100.4°F (38°±/ - 2°C)**. Warm up time: **24 seconds**. Delivery Flow rate: **100mL/min**.

### IV FLUIDS

From a **single battery charge**, capable of warming **1700 mL** of **68°F (20°C) IV fluids** to an output temperature of **100.4°F (38°±/ - 2°C)**. Warm up time: **24 seconds**. Delivery Flow rate: **200mL/min**.

### QUANTUM COMPLETE SYSTEM

ITEM #	NSN #
35-0006	6515-01-693-7810

## TTS-B Thermal Transfusion Set

- Compatible with whole blood, component therapy (including FDP), colloids, and isotonics
- Similar size and weight to traditional IV tubing
- Intelligent sensing & heating elements integrated into the walls of the tubing set
- No disruption to the fluid path or laminar flow of the fluid
- Easy bi-directional connection to the Quantum Controller
- 20 gtt/mL drip rate

### SPECIFICATIONS

TTS-B package (each): 6.5 in. x 11 in.  
Tubing Length: 80 in.  
TTS-B Wgt (each): 2.6 oz  
Packaged Wgt (2 TTS-B & box): 11 oz



TTS-B (2 Pack)

ITEM #	NSN #
35-0005	6515-01-695-4444

## TIS Thermal Infusion Set

- Similar size and weight to traditional IV tubing
- Intelligent sensing and heating elements integrated into the walls of the tubing set
- No disruption to the fluid path or laminar flow of the fluid
- Easy bi-directional connection to the Quantum Controller
- Compatible with standard IV fluids such as colloids and isotonics
- 20 gtt/mL drip rate

### SPECIFICATIONS

TIS packaged (each): 6.5 in. x 11 in.  
Tubing Length: 80 in.  
Weight: TIS (each): 2.1 oz  
Packaged (2 TIS & box): 10 oz

TIS (2 Pack)

ITEM #
35-0004

REV062121

## Crēdo PROMED™ SERIES FOUR CONTAINER

- Temperature-controlled portable medical transport bags that thermally protect the integrity of valuable medical supply payloads for 48 - 96 hours
- Outer Container: Ballistic nylon fabric with durable web strapping and secure buckles
- Featuring the patented Thermal Insulation Chamber (TIC™) panels with phase change material and VIP components contained within to consistently protect medical materials such as blood and platelets within 2-8°C.
- Available in 3 capacities: 2L, 4L & 8L

### CREDO PROMED™ SERIES FOUR CONTAINER

CAPACITY	ITEM #
2 units	35-0011
4 units	35-0012
8 units	35-0013



### SPECIFICATIONS

**2 Units:** Volume: 2L  
Duration: 48 hours  
L 8 in. x W 9 in. x D 10 in.  
Weight: 764 lb

**4 Units:** Volume: 4L  
Duration: 96 hours  
L 10 in. x W 10.5 in. x D 10.5 in.  
Weight: 11.45 lb

**8 Units:** Volume: 8L  
Duration: 72 hours  
L 11 in. x W 15.5 in. x D 11 in.  
Weight: 19.6 lb



*Ideal for use with the Quantum and LifeFlow devices to provide a complete blood delivery package*

### > RESUSCITATE

**500 ml in < 2 minutes  
Stay with Your Patient**

Deliver first bolus quickly to reverse hypotension in patients with shock

- Provide active resuscitation
- Maintain control over delivery

### > RESTORE

**Improve Perfusion /  
Reverse Shock**

Assess using vitals & bedside markers of perfusion, decide if additional fluid bolus is needed.

Early fluids more rapidly reverse shock & **improve outcomes**.<sup>1-3</sup>

### > REASSESS

**Rapid Decision / Disposition**


Each bolus delivered provides an opportunity to evaluate and determine next steps.

Faster resuscitation may lead to:

- Less mechanical ventilation<sup>4</sup>
- More rapid disposition
- **Less donning/doffing of PPE**

## LifeFlow® PLUS BLOOD & FLUID INFUSER

- Hand-operated, rapid infusion device that allows health care providers to deliver fluid & blood products quickly and effectively
- Deliver 500ml of fluid in less than 2 minutes through a 20G IV catheter
- Deliver 1 unit of whole blood in 2 minutes
- Goes from the field to the hospital
- Can be operated with one hand
- Handle and tubing packaged together for quick and easy set-up
- Designed to protect the syringe from bacterial contamination
- Average time to set up and prime: less than 3 minutes



### SPECIFICATIONS

12 in. x 7.5 in. x 2.75 in.  
Weight: 0.95 lb

### LIFEFLOW PLUS BLOOD & FLUID INFUSER

#### ITEM #

30-0287

<sup>1</sup> Lee, S. J., et al. (2014). "Increased fluid administration in the first three hours of sepsis resuscitation is associated with reduced mortality: a retrospective cohort study." *Chest* 146(4): 908-915.

<sup>2</sup> Leisman, D. B., Wie, M., Doerfler, A., Bianculli, M. F., Ward, M., Akerman, J. K., D'Angelo and J. A. Zimmell D'Amore (2016). "Association of Fluid Resuscitation Initiation Within 30 Minutes of Severe Sepsis and Septic Shock Recognition With Reduced Mortality and Length of Stay." *Ann Emerg Med* 68(3): 298-311.

<sup>3</sup> Williams, J. M., J. H. Greenslade, C. A. Dymond, K. Chu, A. F. T. Brown and J. Lipman (2018). "Characteristics, treatment and outcomes for all emergency department patients fulfilling criteria for septic shock: a prospective observational study." *Eur J Emerg Med* 25(2): 97-104.

<sup>4</sup> Leisman D. et al. (2017). "Patterns and Outcomes Associated With Timeliness of Initial Crystalloid Resuscitations in a Prospective Sepsis and Septic Shock Cohort." *Crit Care Med* Oct;45(10):1596-1606.

# Hypothermia Prevention & Management

## HPMK® Hypothermia Prevention & Management Kit

- Includes reinforced Heat Reflective Shell (HRS™) of 4-ply, composite fabric with a protected non-conductive reflective layer
- Built-in hood, integrated fluid absorption pad & top-to-bottom taper providing excellent thermal insulation
- Continuous 1.5 in. Velcro® closures provide rapid 360° access to the casualty with minimal exposure to the elements
- Self-heating, 4-cell liner (Ready-Heat™ Blanket) provides 10 hours of continuous dry heat up to 125°F (52°C)

### SPECIFICATIONS

Shell Packaged:  
H 6.75 in. x W 10.5 in. x D 5.5 in.  
Shell Open: L 78 in. x W 43 in.  
Weight: 3 lb 8 oz

US Patent: 7,766,950 B2



### HPMK - HYPOTHERMIA PREVENTION & MANAGEMENT KIT

ITEM #	NSN #
80-0027	6515-01-532-8056

### HPMK® CONTENTS:



QTY  
**1** HRS®  
Shell

QTY  
**1** Shell  
Liner



NOTE: Litter is sold separately



## Individual First Aid Kits IFAK

### RAPID RESPONSE KIT RTF Edition

- Equipped to provide initial triage of a mass casualty incident and provide treatment for 4 or more patients
- Includes 4 C-A-T® tourniquets with holders and nylon dispenser with 4 rolls of color-coded triage tape

Triage Tape, Set of 4 (Item# ZZ-0591) also available

### RAPID RESPONSE -RTF

COLOR	ITEM #
● Black	85-0716
● Red	85-0717

### SPECIFICATIONS

H 12 in. x W 12 in. x D 6 in.  
Shoulder strap: 28 in. to 50 in.  
Weight: 6 lb 4.5 oz

### RAPID RESPONSE KIT - RTF EDITION:



Red Rapid Response RTF Kit comes standard with orange C-A-T® Tourniquets, red holders and a red tape dispenser



**F**irst responders have increasingly become **the target of violence** in recent history.

Nationwide, **violent crime is rising** at an alarming rate, translating to a higher call volume. With the **sharp decline** in the **police presence** which has historically provided scene safety for Fire/EMS, First Responders are experiencing a **dramatic increase** in **dangerous** and **potentially violent situations** in their everyday operations. Unfortunately rescuers in the Fire/EMS arena can be perceived as just another target for potential life threatening violence along with their law enforcement counterparts.

Now more than ever, your **personal safety is paramount** in the performance of your duties. Wearing **ballistic protective equipment (BPE)** is one of the best ways to **help protect yourself** in a violent environment.

# PROTECT Yourself & SAVE Lives...

Law enforcement transitioned to wearing body armor on a daily basis in the 1980s because they realized it was an essential component of their daily equipment.

**It saves lives!** In the same way, BPE is critical for Rescue Task Force (RTF) and Tactical Emergency Medical Service (TEMS). Fire/EMS should now consider Ballistic Protection Equipment for daily wear on any call that could result in violence or an unsafe environment.

**Protecting Yourself means you can Save others AND Live to serve another day!**



◀ NIJ Level III+ Enhanced Ballistic Rifle Plates



NIJ ▶ 0101.06 Level IIIA Responder PPE Vest Ballistic Protection

▲ Abdominal Extender & Custom ID Patches (each sold separately)



▲ Chest Pouch Medical Kit



Rescue Task Force Rapid Response Kit



▲ MCI WALK



▲ MCI Rapid Response Kit



Explore all your options at **NORTH AMERICAN RESCUE®**


[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

## Customized Trauma Kits

# Customizable Solutions built to meet YOUR requirements

If you haven't found exactly what you need, our team of **subject matter experts** at North American Rescue are standing by to guide you to the **ideal lifesaving solution** tailored to meet all your specifications.

Whether it's the essential gear needed for individuals or customized kits for multiple casualties, a **North American Rescue Specialty Kit** is your answer to the precise aid kit you and your organization need to **increase survivability** should injury occur.



High Quality Tactical Nylon

-OR-

Vacuum Bag



**Call or Email us TODAY to discover all your options for the tailored lifesaving solution YOU can live with!**



Explore all the possibilities at  
**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277





# Individual First Aid Kits IFAK

## MERET EFAK

### MERET EFAK

COLOR	LEVEL	ITEM #
● Black	Basic	80-1200
● Blue	Basic	80-1201
● Red	Basic	80-1202
● Yellow	Basic	80-1203
● Black	Combat Gauze LE	80-1204
● Blue	Combat Gauze LE	80-1205
● Red	Combat Gauze LE	80-1206
● Yellow	Combat Gauze LE	80-1207

- Very compact, hands free medical kit with ample capacity and customizable configuration
- Delivers the tools needed in critical point-of-wounding care for penetrating and other traumatic injuries
- External elastic sleeves secures tourniquet for easy access & rapid deployment
- Provides easy access to gloves and other supplies without opening the main compartment
- Constructed from waterproof, durable, wipe-clean tarpaulin

### MERET EFAK



### SPECIFICATIONS

H 4.5 in. x W 6.25 in. x D 1 in.  
Weight: 0.95 lb



## MERET TFAK

### MERET TFAK

COLOR	LEVEL	ITEM #
● Black	Basic	80-1208
● Blue	Basic	80-1209
● Black	Combat Gauze LE	80-1210
● Blue	Combat Gauze LE	80-1211



Same features as the MERET EFAK but in a larger bag with an additional C-A-T Tourniquet and extra pair of Responder gloves

### MERET TFAK



### SPECIFICATIONS

H 4.5 in. x W 8 in. x D 1.5 in. | Weight: 1.26 lb

# Multiple Injury Trauma Kits

## Mass Casualty Incident **MCI WALK**



- All-in-one trauma treatment & evacuation platform for treatment of multiple trauma patients
- Choose from two different litter configurations: 8 non-rigid, disposable QuikLitters or 1 full-size, quad-fold, lightweight MedEVAC Litter
- Pocket system holds self- or buddy-aid point-of-wounding medical supplies
- Bag features back pack straps for easy transport from the vehicle to the point of need
- Ideal for use as a vehicle trauma kit

(80-0565/80-0458) Red MCI WALKs contain orange C-A-T® Tourniquets

(#80-0458) Red MCI WALK with MedEVAC Litter contains an orange MedEVAC Litter

### MCI WALK

BAG COLOR	LITTER OPTION	C-A-T	ITEM #
● Black	8 QuikLitters	Black	80-0564
● Red	8 QuikLitters	Orange	80-0565
● Black	MedEVAC Litter (Black)	Black	80-0457
● Red	MedEVAC Litter (Orange)	Orange	80-0458

### SPECIFICATIONS

H 22.75 in. x W 16.5 in. x D 11.5 in.  
 Weight with 8 QuikLitters: 18 lb 12 oz  
 Weight with MedEVAC Litter: 29 lb 15 oz

### MCI WALK with QUIKLITTERS:

 QTY 8 C-A-T® Tourniquet	 QTY 8 ETD™ 6 in. Dressing	 QTY 5 Combat Gauze® LE	 QTY 8 Gauze, Compressed	 QTY 5 28F Naso Airway & Lube
 QTY 8 HyFin® Vent. Twin Pk	 QTY 8 NAR QuikLitter	 PAIR 8 Nitrile Gloves, Lg	 QTY 20 Combat Casualty Card	 QTY 5 Permanent Marker, Large

### MCI WALK with MEDEVAC LITTER:

 QTY 8 C-A-T® Tourniquet	 QTY 8 ETD™ 6 in. Dressing	 QTY 5 Combat Gauze® LE	 QTY 8 Gauze, Compressed	 QTY 5 28F Naso Airway & Lube
 QTY 8 HyFin® Vent. Twin Pk	 QTY 1 Rigid Litter MedEVAC	 PAIR 8 Nitrile Gloves, Lg	 QTY 20 Combat Casualty Card	 QTY 5 Permanent Marker, Large



Bag only option also available  
at NARescue.com



online video  
available

< scan to view

## Splinting & Patient Evacuation

# SLISHMAN TRACTION SPLINT



STS-  
STANDARD

STS-  
COMPACT

Trusted by major EMS agencies including Boston EMS, Dallas Fire-Rescue, LA County Fire, and FDNY

- One of the smallest and lightest traction splints on the market
- One size fits all, from small children to large adults
- Can be applied in under a minute with minimal training

- Single connected device—no loose parts, which may get lost
- Application possible even in cases of lower leg/foot amputation
- Radiolucent

### SLISHMAN TRACTION SPLINT

OPTION	ITEM #
STS	50-0041
STS-Compact	50-1008

### SPECIFICATIONS

**Standard:** Packaged: L 23 in. x W 3 in. x D 3 in.  
Fully Extended: L 45.5 in. x Diameter 0.75 in.  
Weight: 21 oz

**Compact:** Packaged: L 13 in. x W 3 in. x D 3 in.  
Fully Extended: L 45.5 in. x Diameter 0.75 in.  
Weight: 13 oz

## PELVIC BINDER from CRO Medical

- BOA Precision Dial System
- Hook and loop retention for Pressure Delivery Device
- Approximately half the size in cubic inches of other binders
- Lightweight 4 in. double-back wrap
- Easy, one-hand tightening
- Easily achieves at least 33 lb (150N) of circumferential force
- Easily reduces and stabilizes open-book pelvic fractures using its patented tightening system



### PELVIC BINDER

ITEM #

50-0044

### SPECIFICATIONS

Standard Waist Size: 32 in. - 50 in. hips  
Weight: 0.75 lb

## SAM® PELVIC SLING for Pelvic Ring Fractures

- Force-controlled, ideal compression every time without the guesswork
- Rapid, simple, three-step application, usually in less than one minute



### SPECIFICATIONS

H 7.25 in. x W 7.75 in. x D 4 in.  
Hip Circumference: 32 in. to 50 in.  
Weight: 9.8 oz

### SAM PELVIC SLING

ITEM #	NSN #
50-0004	6515-01-509-6866

## NAR QUIKLITTER LITE

- Emergency carry disposable litter that is designed for rapid deployment for the movement of casualties for mass casualty incidents, disasters, active shooter events, as well as for everyday patient movements
- Made of lightweight, heavy-duty nylon
- Features ten large carry handles constructed of heavy-duty nylon webbing
- Load capacity tested to 750 lb
- Vacuum-packed to conserve space and fit in a cargo pocket



### SPECIFICATIONS

Payload 750 lb

Litter Bed:  
L 74 in. x W 27 in.

Weight: 12 oz

### QUIKLITTER LITE

ITEM #

60-0071

## UNO EMERGENCY SPLINT

- Modular Design
- Stepless Adjustment
- Small and Lightweight
- Compatible with TCCC Equipment
- Pack of 3 includes additional ankle strap for stabilizing the ankle and enables the device to provide traction



### SPECIFICATIONS

Packaged - L 21.50 in. x W 10.55 in. x D 1.50 in.  
Unpackaged - L 15.75 in. x W 2.17 in. x D 2.36 in.  
Weight: 0.91 lbs

### UNO EMERGENCY SPLINT

ITEM #	OPTION
50-0045	Single
50-0046	Pack of 3

Take your training to the  
**NEXT LEVEL**  
with superior **REALISM**

Tactical  
Operation  
Medical Manikin  
**SERIES**

**TOMM**Manikin™

**BASIC | BLAST | GSW | WATER | CONFINED SPACE | BURN | CBRNE | FEMALE**

**Interactive Animatronics:**

- Lifelike Breathing Mechanisms
- Functional Anatomical Airway
- Variable Radial Pulse Rates & Pressure
- High Fidelity Bleeding from Multiple Wounds
- Customizable, Multi-language Voice Tracks
- Rugged Softech® Skin in a variety of skin tones

**Wireless Controls:**

- MILSPEC Wireless Remote Instructor Tablet & Software
- Stand alone point-to-point communications with .5 mile line of site range

**Flexible Wound Configurations:**

**Basic | Blast | GSW | Water | Confined Space | Burn | CBRNE**  
Custom configurations including Female Anatomy on request

**Medical Procedural Capabilities:**

**MASSIVE HEMORRHAGE:**

- Tourniquet Application (proximal, junctional & abdominal)
- Active Bleeding Control by Direct Pressure on Pressure Points
- Wound Packing

**ADVANCED AIRWAY:**

- Nasopharyngeal/Oral/Supraglottic Airway Insertion
- Endotracheal Intubation
- Surgical Airway and Tissue Dissection with Realistic Bleeding

**RESPIRATIONS:**

- Chest Seal Application
- Bilateral Needle Decompression
- Chest Tube Capable

**CIRCULATION:**

- IO (sternal & humoral)
- Antecubital IV Access on uninjured upper limbs

REV050121



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

For more info, contact **NAR Simulation Sales** | [manikins@narescue.com](mailto:manikins@narescue.com) | **866.839.3455**

## ARMADILLO Medication & Storage Case

- Compact solution for keeping critical medications safe, secure and easily accessible
- Features Automatic Pressure Equalization Valve
- Hardware constructed of Marine grade 316 SS
- Case is water-resistant and crush-proof
- Hook & loop attachment system allows storage compartment to be tailored to hundreds of custom configurations

### SPECIFICATIONS

Ext: L 7.95 in. x W 5.1 in. x D 2.14 in.  
Int: L 6.45 in. x W 3.87 in. x D 1.75 in.  
Weight: 10.5 oz

### ARMADILLO MEDICATION CASE \* NATO NSN

LID STYLE	ITEM #	NSN #
● Black	ZZ-0081	8145-01-573-2533
○ Clear	ZZ-0259	6545-01-594-9399*



Clear Lid  
ZZ-0259 only  
includes Loop Liner  
on the bottom

**also available:** Half Inch Hook Discs with pressure sensitive adhesive  
Set of 250 • Item# ZZ-0258 • NSN# 8145-01-617-9075\*

## SCISSOR LEASH™

- Salt water resistant
- Self cleaning
- Helps keep track of your shears

### SPECS

H 6 in. x W 16 in. x D 0.6 in.  
Cable Length - 36 in.  
Weight: 1.6 oz

### SCISSOR LEASH

ITEM #	NSN #
70-0005	6515-01-540-7226



## 6 oz MOLLE Retractable GEAR LEASH

- MOLLE Mount compatible with 360° rotation
- High break strength / Salt water resistant
- 36 in. Extension / Spectra Nylon Line
- Self-Flushing
- Stainless Steel Spring & Hardware

### GEAR LEASH

ITEM #
70-0155



### SPECS

Retracted:  
H 3.5 in. x W 1.6 in. x D 1.0 in.  
Cable Length - 36 in.  
Weight: 2.0 oz

## Blue RESPONDER NITRILE GLOVES

- Form-fitting, ambidextrous design with textured fingertips
- Thick palm ensures durability
- Textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves
- Available in pre-packaged pairs (kits) or in bulk (box)

COLOR:  
Blue

### SPECIFICATIONS

Bag (25 pair): L 4.75 in. x W 6 in. x D 4 in.  
Weight: 1 lb 1 oz

Box (Bulk): L 11 in. x W 5 in. x D 3.75 in.  
Weight: 2 lb 4 oz



### Blue Responder BULK

SIZE	PKG   QTY	ITEM #	SIZE	PKG   QTY	ITEM #
S	Box   50 pair	70-0315	L	Box   50 pair	70-0317
M	Box   50 pair	70-0316	XL	Box   50 pair	70-0318

### Blue Responder PAIRS

SIZE	PKG   QTY	ITEM #
S	Bag   25 pair	ZZ-1059
M	Bag   25 pair	ZZ-1060
L	Bag   25 pair	ZZ-1043
XL	Bag   25 pair	ZZ-1061

\* also available in black at [NARescue.com](http://NARescue.com)

Heavy Duty Trauma **XSHEAR®**

- Hardened stainless steel blades that are twice as thick as most trauma shears, sharpened to a razor edge
- Patented blunt tip and curved design for gentler edge near the skin
- Soft touch inner handles for comfort and slip resistance
- Heavy duty center bolt that won't loosen over time

**XSHEAR**

HANDLE	BLADE	ITEM #
● Black	Black	70-0290
● Blue	Silver	70-0291
● Orange	Black	70-0292
● Pink	Black	70-0293
● Red	Black	70-0294



**SPECIFICATIONS**

L7.5 in x 3.5 in. x 0.118 in.  
Cutting area length: 2 in.  
Weight: 4 oz



Shears & Holster are sold separately

**RIGID XSHEAR® HOLSTER**

**ITEM #**

70-0295

**Rigid XSHEAR® HOLSTER**

- Custom designed for the XSHEAR® Heavy Duty Trauma Shears
- MOLLE Compatible adjustable clip
- Flex lock design securely holds shears in place with a positive click
- Universal design allows shears to be inserted either direction
- Rugged injection molded ABS Plastic

**SPECIFICATIONS**

L 5 in. x W 3.75 in. x D 1.75 in. | Weight: 3.2 oz

**Soft XSHEAR® HOLSTER**

- Designed for XSshears but compatible with most shears
- Durable Spandex over soft foam with nylon interior
- Powder coated heavy steel belt clip
- Allows Left or Right hand carriage



**SOFT XSHEAR® HOLSTER**

**ITEM #**

70-0314

**SPCS**

L 7 in. x W 3.5 in. x D 0.5 in.  
Weight: 3.2 oz



**SPECIFICATIONS**

Shear Length: 19 in.  
Closed Length: 5 in.  
Weight: 5.8 oz

**Leatherman® RAPTOR™ MULTI-TOOL**

- Integrated multi-tool specific to emergency medical situations
- High quality 420HC stainless steel medical shears
- Foldable shears come apart for easy cleaning
- Injection-molded polymer Holster on a rotating clip; MOLLE-compatible

**RAPTOR MULTI-TOOL**

HANDLE	ITEM #	NSN #
● Black	70-0014	5110-01-626-8766
● Blue	70-0264	----
● Coyote	70-0265	----
● Red	70-0266	----
● Green	70-0326	----
● Org/Blk	70-0327	----



# NORTH AMERICAN RESCUE®

Education and Training

## North American Rescue Education and Training

has roots that run deep and wide in the world of pre-hospital medicine and technical rescue. Our team was built from the ground up on the vision and experience of seasoned military medics seeking to fill the training gaps commonly encountered in the operational environment.

Since then, we have continued to broaden our vision to include other invaluable first responder perspectives such as pre-hospital healthcare providers in civilian, law enforcement, emergency medical services and community public safety initiatives.


## CORE Competencies

- Continuing Medical Education (CMEs)
- Prolonged Casualty Care
- Austere Medicine
- Technical Rescue
- Emergency Medical Services
- Critical Care Flight Medicine



## NAR Training

delivers a broad range of curricula which frequently include **Hands-on Skills Stations** and **Field Scenarios** designed to develop both competence and confidence. Using **Professional Role Players, Moulage Artists, Subject Matter Experts, High Fidelity Manikins, Re-animated Cadavers** and **State of the Art Facilities**, our team delivers unparalleled medical instruction to today's heroes.



## Training Locations:

- NAR Training
- Training Partners & Affiliates

## ACCREDITATIONS:



American Heart Association®

Now SAME DAY Cards Issued At End Of Course!



Explore training opportunities online

**NAReducation.com** ➤ Phone: 888.689.6277 ➤ Email: [NARTraining@NARescue.com](mailto:NARTraining@NARescue.com)

888-689-6277 • [www.NARescue.com](http://www.NARescue.com) | 23

MASSIVE HEMORRHAGE | AIRWAY | RESPIRATION | CIRCULATION | HYPOTHERMIA



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

REV092022



*Explore all your options online...*