

TOURNIQUETS | FIRST AID KITS | TRAUMA KITS | RESCUE TOOLS | LITTERS

LE

LAW ENFORCEMENT
Catalog

**CALL FOR
DISCOUNTED PRICING**

888.689.6277

Shop ALL products at www.NARescue.com



Since 1996, North American Rescue has focused on **increasing survivability** on both the battlefield and the homefront. With our innovative medical/rescue equipment, NAR continues to lead the way by delivering the **tailored, mission-critical healthcare solutions** needed by today's first responders to counter the rapidly evolving threats of the 21st century. See what a lifesaving difference our combat proven make in the prehospital environment!

CONTENTS

TOURNIQUETS		3-4
FIRST AID KITS		4 -12
TRAUMA KITS		13-16
EXTRAS		17



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277





CUSTOMIZABLE SOLUTIONS

BUILT TO MEET YOUR REQUIREMENTS

If you haven't found exactly what you need, our team of *subject matter experts* at North American Rescue are standing by to guide you to the *ideal lifesaving solution* tailored to meet all your specifications.

Whether it's the essential gear needed for individuals or customized kits for multiple casualties, a *North American Rescue Specialty Kit* is your answer to the precise aid kit you and your organization need to *increase survivability* should injury occur.



Call or Email us TODAY to discover all your options for the tailored lifesaving solution YOU can live with!



Explore all the possibilities at
NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277



Massive Hemorrhage | Tourniquets

C-A-T® Combat Application Tourniquet Gen 7

- **Single route buckle** decreases application time to help decrease blood loss, allows for slack removal for proper TQ application
- **Applied same way** to both upper and lower extremities
- True **one-handed tourniquet**
- **Scientifically proven** 100% effective in occluding blood flow in both upper and lower extremities by the USAISR and recommended by CoTCCC

The FASTEST, SAFEST, Most EFFECTIVE Prehospital Field Tourniquet!



SPECIFICATIONS

Packaged: H 6.5 in. x W 2.4 in. x D 1.5 in.
Device: L 37.5 in.
Weight: 2.7 oz

COMBAT APPLICATION TOURNIQUET

COLOR	STYLE	ITEM #	NSN #
Black	Tactical	30-0001	6515-01-521-7976
Orange	Rescue	30-0023	6515-01-612-7129
Blue	Trainer	30-0033	6910-01-560-2972



online video available

< scan to view

RIGID C-A-T® TQ CASE with Shirt Shield

- Allows you to carry a Gen 7 C-A-T® Tourniquet on your duty, patrol, or pistol belt so you have life-saving equipment on your person at all times
- Includes a Shirt Shield complete with mounting hardware to protect your uniform shirt from undue wear caused by the friction from the windlass clip
- Blade-Tech TekLok Belt Attachment fits up to a 2.25 in. belt
- Can be mounted either vertically or horizontally or attached to a Safariland Drop Leg Shroud
- Injection Molded Nylon Case stands up to rugged LE/MIL usage
- Cross medical identifier patch included

RIGID C-A-T TQ CASE with Shirt Shield

COLOR	ITEM #
Black	80-0944
Black, Basket Weave	80-0962
Black, Polished	80-0963
Coyote	80-0945
Ranger Green	80-0946
Navy	80-1161

*More colors available at NARescue.com

SPECIFICATIONS

H 5.88 in. x W 2.5 in. x D 2 in. | Weight: 5.36 oz



NOTE: All Tourniquets sold separately

Nylon C-A-T® HOLDER

- Tailored specifically for your C-A-T®, to protect it from exposure to the elements
- PALS/MOLLE-style connectors for versatile attachment to the duty vest, gear or belt (mounted horizontally or vertically) for rapid & easy access
- Easy-open pull elastic tab with tourniquet ID patch – “TQ”



NYLON TQ HOLDER

COLOR	ITEM #
● Black	30-0057
● Coyote	30-0044
● Multi-Cam	30-0059
● Red	30-0128
● Ranger Green	30-0208
● RealTree®	30-0207
● OD Green	30-0053

SPECIFICATIONS

L 6.25 in. x W 1.75 in. x D 1.5 in.
Weight: 2.17 oz

LASER CUT C-A-T HOLDER

- Constructed with laser cut laminated Cordura® nylon
- Easy Pull Elastic Cover
- PALS/MOLLE-style connections
- Accommodates GEN 6 and GEN 7 C-A-T® Tourniquets
- Weighs less than traditional TQ holders

SPECIFICATIONS

L 7.25 in. x W 2.75 in. x D 1.5 in.
Weight: 1.29 oz

LASER CUT TQ HOLDER

COLOR	ITEM #
● Black	30-0297
● Coyote	30-0298
● Ranger Green	30-0299
● OCP	30-0306



Individual First Aid Kits IFAK

EDC *Every Day Carry* ANKLE TRAUMA KIT

- Very lightweight, low profile, versatile carrier that allows for carry of essential trauma and EDC equipment with quick, single-handed, ambidextrous access to components
- Pocket designed to accommodate most firearms magazines including double stack magazines
- Separate concealed compartment for essential EDC items or ID
- Designed to accommodate a C-A-T® Tourniquet and gauze/hemostatic dressing/bandage
- Advanced materials that provide user comfort and long-term wearability
- Individual pouches for each carried component with strong hook and loop closures
- One size fits most with convenient hook and loop adjustability
- Available with & without medical equipment

EDC ANKLE TRAUMA KIT

OPTIONS	ITEM #
Wound Packing	80-1157
BCD	80-1158
Carrier Only	30-0294T

SPECIFICATIONS

L 17.75 in. x H 4 in. x D 0.25 in.
Weight: 0.657 lb



EVERY DAY CARRY ANKLE TRAUMA KIT

 QTY 1 C-A-T® Tourniquet	 QTY 1 HyFin® Vent Compact, Twin Pk	 PAIR 1 Responder Nitrile Gloves, Lg	 QTY 1 Gauze Wound Packing	 QTY 1 Bleeding Control Dressing
---	--	---	---	---

IFAK Individual First Aid Kits



Individual Patrol Officer Kit **IPOK**


- Individual hemorrhage control kit
- Compact, vacuum-sealed & packaged to fit into a pocket, vest pouch or individual bag
- Rugged, durable pouch features NAR's signature Red Tip Technology® ensuring quality & performance

SPECIFICATIONS

H 4 in. x W 6.5 in. x D 2.75 in.
Weight: 7.4 oz



IPOK:



QTY **1** C-A-T® Tourniquet

QTY **1** ETD™ 4 in. Flat

QTY **1** Choice of Wound Pk Gauze or Hemostatic Dressing

PAIR **1** Nitrile Gloves, Lg

IPOK - INDIVIDUAL PATROL OFFICER KIT

HEMOSTATIC/GAUZE	ITEM#
Wound Packing Gauze	80-0167
Combat Gauze, LE	80-0168

ETAK

Enhanced Trauma Aid Kit

Designed to pair with a separate C-A-T® Tourniquet & holder (not included) for a great option to add more capability to previously owned kit

- Compact IFAK platform accessible from a belt, duty belt, assault vest/body armor or any equipment with PALS/MOLLE attachment points.
- Featuring multiple interior elastic loops keep your kit organized and secure
- Exterior tear away panel is easily removable for added mobility and access to your kit
- MED patch included

Every ETAK includes these core items:



QTY **1** ETD™ 4 in. Mini Responder

QTY **1** HyFin® Vent Compact, Twin

PAIR **1** Nitrile Gloves, Lg



Bag only option also available at NARescue.com

PLUS+

BASIC also includes:



QTY **1** Wound Packing Gauze


QTY **1** Permanent Marker, Small

QTY **1** Combat Gauze® LE

QTY **1** Permanent Marker, Small

with **COMBAT GAUZE** also includes:

ADVANCED also includes:



QTY **1** Combat Gauze® LE

QTY **1** 28F Naso Airway & Lube

QTY **1** ARS® Needle, 10g



SPECIFICATIONS

H 5.75 in. x W 4 in. x D 1.75 in.
Weight: 0.65 lb

ETAK ENHANCED TRAUMA AID KIT

COLOR	LEVEL	ITEM #
● Black	Basic	80-1168
● Coyote	Basic	80-1169
● Ranger Green	Basic	80-1171
● MultiCam	Basic	80-1170
● Black	Combat Gauze LE	80-1172
● Coyote	Combat Gauze LE	80-1173
● Ranger Green	Combat Gauze LE	80-1175
● MultiCam	Combat Gauze LE	80-1174
● Black	Advanced	80-1176
● Coyote	Advanced	80-1177
● Ranger Green	Advanced	80-1179
● MultiCam	Advanced	80-1178

Individual First Aid Kits IFAK

IPOK HOLDER

- Protects your medical kit from the elements
- Enclosure can be secured by Fastex snap buckle, hook-and-loop or both



IPOK HOLDER

COLOR	ITEM#
● Black	80-0205T
● Coyote	80-0207
● OD Green	80-0254T
● Red	80-0358T

SPECIFICATIONS

H 6.5 in. x W 6 in. x D 3.5 in.
Weight: 4.1 oz

VISOR MED POUCH

- Securely mounts to a vehicle visor
- Orange webbing handle to grab for easy deployment
- Pouch with easy-open buckle



VISOR MED POUCH

COLOR	ITEM#
● Black	80-0269
● Coyote	80-0270
● Red	80-0360

SPECIFICATIONS

H 6.5 in. x W 5.5 in. x D 3 in.
Weight: 5.4 oz

HEADREST IFAK

- Designed to provide personnel with an easily accessible, quick to deploy medical pouch that conveniently mounts to a vehicle's headrest
- Fully loaded with emergency medical equipment for point-of-wounding treatment until definitive care is available
- Durable nylon construction features multiple interior elastic loops allowing for the contents to be organized and secured
- Easy to grab side handles allow for easy deployment
- IFAK pouch only is also available
- The pouch design allows for the contents of an individual first aid kit or vacuum sealed kit

SPECIFICATIONS

H 5.25 in. x
W 7.25 in. x D 4.25 in.
Weight: 1.63 lb



HEADREST IFAK



HEADREST IFAK

COLOR	OPTIONS	ITEM #
● Black	Basic	80-1091
● Coyote	Basic	80-1150
● Black	Combat Gauze LE	80-1151
● Coyote	Combat Gauze LE	80-1152
● Black	Bag Only	80-1090T
● Coyote	Bag Only	80-1149T

IFAK Individual First Aid Kits

M-FAK® Mini First Aid Kit - LE

- Contains the essential point-of-wounding equipment needed for treating penetrating, blast and other traumatic injuries
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access

Bag only option also available at NARescue.com

SPECS

H 6 in. x W 3 in. x D 3 in.
Weight: 13 oz



LE M-FAK | ROO M-FAK BASIC or BASIC with Combat Gauze:

BASIC	-OR-	BASIC with Combat Gauze	-AND-				
QTY 1 Gauze S-Rolled		QTY 1 Combat Gauze® LE	QTY 1 C-A-T® Tourniquet	QTY 1 ETD™ 4 in. Flat	QTY 1 HyFin® Vent Compact, Twin	PAIR 1 Nitrile Gloves, Lg	

LE M-FAK | ROO M-FAK ADVANCED:

QTY 1 C-A-T® Tourniquet	QTY 1 ETD™ 4 in. Flat	QTY 1 Combat Gauze® LE	QTY 1 28F Naso Airway & Lube	QTY 1 HyFin® Vent Compact, Twin	QTY 1 ARS® Needle, 14g	PAIR 1 Nitrile Gloves, Lg
-------------------------	-----------------------	------------------------	------------------------------	---------------------------------	------------------------	---------------------------

M-FAK MINI FIRST AID KIT - LE

COLOR	LEVEL	ITEM #
Black	Basic	80-0494
Coyote	Basic	80-0495
Gray	Basic	80-0872
Multi-Cam	Basic	80-0496
OD Green	Basic	80-0497
Ranger Green	Basic	80-0996
Black	Combat Gauze LE	80-0498
Coyote	Combat Gauze LE	80-0499
Gray	Combat Gauze LE	80-0873
Multi-Cam	Combat Gauze LE	80-0500
OD Green	Combat Gauze LE	80-0501
Ranger Green	Combat Gauze LE	80-0998
Black	Advanced	80-0502
Coyote	Advanced	80-0503
Gray	Advanced	80-0874
Multi-Cam	Advanced	80-0504
OD Green	Advanced	80-0505
Ranger Green	Advanced	80-0997



SPECS

H 7 in. x W 3.5 in. x D 2.25 in.
Weight: 6.4 oz

Compact Officer Response Emergency Kit CORE KIT

Designed to go where you go, this compact medical kit will fit in a BDU pocket, vest pouch or Go-Bag and provides personnel with essential life-saving bleeding control equipment that can be kept on their person. The CORE Kit is packaged in a compact, rugged, clear vacuum-sealed bag with easy-to-open tear notches.

- Compact design with the essential life-saving equipment needed to address moderate and severe traumatic bleeding
- Contains the #1 CoTCCC Recommended Prehospital Tourniquet and Hemostatic Gauze
- The Best Tourniquet paired with the Best Hemostatic Gauze and an effective pressure bandage at the Best Price
- Meets the Tactical Emergency Casualty Care (TECC) Guidelines for hemorrhage control

CORE KIT:

QTY 1 C-A-T® Tourniquet	QTY 1 ETD™ 4 in. Responder, Flat	QTY 1 Combat Gauze® LE
-------------------------	----------------------------------	------------------------

CORE KIT

ITEM #
85-0834

ROO M-FAK *Mini First Aid Kit*

Rapid one handed access to your tourniquet without ever having to open the full kit!

Contains the same packout options as the M-FAK: **Basic**, **Basic with Combat Gauze LE** and **Advanced** (see page 28)

ROO M-FAK MINI FIRST AID KIT

COLOR	LEVEL	ITEM #
● Black	Basic	80-1049
● Coyote	Basic	80-1050
● OD Green	Basic	80-1051
● Ranger Green	Basic	80-1052
● Black	Combat Gauze LE	80-1053
● Coyote	Combat Gauze LE	80-1054
● OD Green	Combat Gauze LE	80-1055
● Ranger Green	Combat Gauze LE	80-1056
● Black	Advanced	80-1057
● Coyote	Advanced	80-1058
● OD Green	Advanced	80-1059
● Ranger Green	Advanced	80-1060



L-FAK *Lumbar First Aid Kit*

- Durable 500-denier exterior carrier
- Compact design to reduce space needed for life-saving equipment
- Designed for wear in the small of the back
- Multiple interior elastic loops for securing gear
- Rapid access to medical gear
- C-A-T® tourniquet sold separately



Every L-FAK includes these core items:



PLUS+ BASIC
also includes:



COMBAT GAUZE™
also includes:



ADVANCED also includes:



SPECIFICATIONS

Pkgd - H 4.75 in. x W 5.75 in.
Weight: 0.88 lbs
Adv. Weight: 0.94 lbs

LUMBAR FIRST AID KIT

COLOR	LEVEL	ITEM #
● Black	Basic	85-3627
● Coyote	Basic	85-3633
● Ranger	Basic	85-3630
● Black	w/ Combat Gauze™	85-3628
● Coyote	w/ Combat Gauze™	85-3634
● Ranger	w/ Combat Gauze™	85-3631
● Black	Advanced	85-3629
● Coyote	Advanced	85-3635
● Ranger	Advanced	85-3632

Bag only option also available at NARescue.com

REFLEX SYSTEM ROLL & CARRIER

- Primary construction of heavy-duty nylon laminate allows rapid deployment of medical supplies
- Features High Abrasion-Resistant Neoprene Kevlar® composite grab handles for long-term durability
- Can mount vertically or horizontally on belts 1.5 in.- 2.25 in., and horizontally on MOLLE
- Carry handle and shock cord loops allow the roll to be attached to a patient, or hung from an I.V. rack
- "Mini" MOLLE allows for precise fit on belts as narrow as 1.5 in. or two-row MOLLE panels
- Features loop panels on med roll and carrier for label patches (Medical Cross Patch included)
- Medical contents sold separately



REFLEX SYSTEM ROLL & CARRIER

COLOR	ITEM #
● Black	80-1117
● LE Blue	80-1118
● Coyote	80-1119
● Gray	80-1120
● Multi-Cam Blk	80-1121
● Multi-Cam	80-1122
● OD Green	80-1123

SPECIFICATIONS

L 8 in. x W 3 in. x H 5.5 in.
Weight: .58 lb

Two-piece IFAK system designed to carry organized supplies with ambidextrous accessibility



▲
Carrier

▲
Roll

Reflex System VEHICLE MOUNT



REFLEX MOUNT

COLOR	ITEM #
● Black	ZZ-1368
● Coyote	ZZ-1369
● Gray	ZZ-1370

- Attaches ReFlex™ Med Roll to the back of a vehicle headrest
- Fits virtually all vehicles with detachable headrests
- Elastic cinch strap allows instant retrieval of the Med Roll from either side of the vehicle
- All-nylon construction with built-in elastic tourniquet loop

SPECIFICATIONS

L 5.5 in. x W 0.125 in. x H 6 in.
Weight: 0.21 lb



REFLEX SYSTEM IFAK

COLOR	LEVEL	WGT	ITEM #
● Black	Basic	1.58 lb	80-1127
● LE Blue	Basic	1.58 lb	80-1128
● Coyote	Basic	1.58 lb	80-1129
● Gray	Basic	1.58 lb	80-1130
● Multi-Cam Black	Basic	1.58 lb	80-1131
● Multi-Cam	Basic	1.58 lb	80-1132
● OD Green	Basic	1.58 lb	80-1133
● Black	Combat Gauze	1.6 lb	80-1134
● LE Blue	Combat Gauze	1.6 lb	80-1135
● Coyote	Combat Gauze	1.6 lb	80-1136
● Gray	Combat Gauze	1.6 lb	80-1137
● Multi-Cam Black	Combat Gauze	1.6 lb	80-1138
● Multi-Cam	Combat Gauze	1.6 lb	80-1139
● OD Green	Combat Gauze	1.6 lb	80-1140
● Black	Advanced	1.61 lb	80-1141
● LE Blue	Advanced	1.61 lb	80-1142
● Coyote	Advanced	1.61 lb	80-1143
● Gray	Advanced	1.61 lb	80-1144
● Multi-Cam Black	Advanced	1.61 lb	80-1145
● Multi-Cam	Advanced	1.61 lb	80-1146
● OD Green	Advanced	1.61 lb	80-1147

REFLEX SYSTEM IFAK

- Developed with direct input from active-duty medical personnel, the Reflex System Roll & Carrier fully loaded with essential medical equipment to treat life threatening injuries in the field of operation
- Intentional organization options put all medical supplies within easy reach
- Primary construction of heavy-duty nylon laminate allows rapid deployment of medical supplies
- Features High Abrasion-Resistant Neoprene Kevlar® composite grab handles for long-term durability
- Can mount vertically or horizontally on belts 1.5 in.- 2.25 in., and horizontally on MOLLE
- Carry handle and shock cord loops allow the roll to be attached to a patient, or hung from an I.V. rack
- "Mini" MOLLE allows for precise fit on belts as narrow as 1.5 in. or two-row MOLLE panels



SPECIFICATIONS

L 8 in. x W 3 in. x H 5.5 in. | Weight: 1.6 lb

Every REFLEX SYSTEM IFAK includes these core items: //

 QTY 1 C-A-T® Tourniquet	 QTY 1 ETD™ 4 in. Flat	 QTY 1 Wound Packing Gauze	 QTY 1 HyFin® Vent, Twin Pk	 QTY 1 Survival Blanket
 QTY 1 Surgical Tape 2" Flatfold	 QTY 1 Trauma Shears, Lg	 PAIR 2 Nitrile Gloves, Lg	 QTY 2 Combat Casualty Card	 QTY 1 Permanent Marker, Large

PLUS+ BASIC

also includes:

 QTY 1 Wound Packing Gauze
--

COMBAT GAUZE

also includes:

 QTY 1 Combat Gauze™ Z-fold
--

ADVANCED also includes:

 QTY 1 Combat Gauze™ Z-fold	 QTY 1 28F Naso Airway & Lube	 QTY 1 ARS® Needle, 10g
--	--	--

RIG SERIES EAGLE IFAK®

- Tourniquet attached for quick deployment by banding it to inside of front opening flap
- Protective hood with quick pull latch release helps protect IFAK components and includes sewn in embroidered medical cross identifier
- Detachable (hook & loop backed) shock-cord inner panel adjusts to secure internal components and is easily removable for rapid access

Bag only option
also available at
NARescue.com



SPECIFICATIONS

H 8 in. x W 4 in. x D 3 in.
Kit Packed Weight: 118 lb
Empty Bag Weight: 7 oz

RIG SERIES EAGLE IFAK

COLOR	LEVEL	ITEM #
● Black	Basic	80-0509
● Coyote	Basic	80-0510
● Multi-Cam	Basic	80-0511
● Ranger Green	Basic	80-0991
● Black	Combat Gauze	80-0517
● Coyote	Combat Gauze	80-0518
● Multi-Cam	Combat Gauze	80-0519
● Ranger Green	Combat Gauze	80-0993
● Black	Advanced	80-0513
● Coyote	Advanced	80-0514
● Multi-Cam	Advanced	80-0515
● Ranger Green	Advanced	80-0992

RIG SERIES EAGLE BASIC:



RIG SERIES EAGLE with COMBAT GAUZE LE:



RIG SERIES EAGLE ADVANCED:



Individual First Aid Kits IFAK

T.O.R.K. TACTICAL OPERATOR RESPONSE KIT

COLOR	LEVEL/HEMOSTATIC	ITEM #
Black	Basic	85-0008
Coyote	Basic	85-0009
Digital Camo	Basic	85-2300
Multi-Cam	Basic	85-0012
Ranger Green	Basic	85-2301
Red	Basic	85-1805
Black	Basic w/Combat Gauze	80-1097
Coyote	Basic w/Combat Gauze	80-1098
Digital Camo	Basic w/Combat Gauze	80-1099
Multi-Cam	Basic w/Combat Gauze	80-1100
Ranger Green	Basic w/Combat Gauze	80-1102
Red	Basic w/Combat Gauze	80-1108
Black	Advanced w/Combat Gauze	85-0048
Coyote	Advanced w/Combat Gauze	85-0049
Digital Camo	Advanced w/Combat Gauze	85-0050
Multi-Cam	Advanced w/Combat Gauze	85-0052
Ranger Green	Advanced w/Combat Gauze	85-2302
Red	Advanced w/Combat Gauze	85-1806
7 options online	Advanced w/ChitoGauze	NARescue.com

T.O.R.K. Tactical Operator Response Kit

- Streamlined, compact IFAK designed to MOLLE to the assault vest or body armor
- Includes essential life-saving equipment to be used at the point of wounding
- For self-aid/buddy-aid injuries caused by penetrating or blast trauma

T.O.R.K. ReSUPPLY KITS*

CONTENT OPTIONS	ITEM #
Basic	80-1037
Basic Combat Gauze	80-1103
Advanced Combat Gauze	80-1039
Advanced ChitoGauze	80-1038



*ReSupply Kits include contents only. Nylon bag is sold separately.

T.O.R.K. BASIC/BASIC with Combat Gauze:



SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.75 in.
Weight: 1 lb 5 oz

T.O.R.K. ADVANCED:



Bag only options also available at NARescue.com

Trauma Kits

MINI MEDIC KIT

- Compact, multi-purpose trauma aid bag with vital medical equipment for point of wounding treatment of multiple injuries
- Two external pockets with two external tourniquet holders
- One external pouch with additional drop-down pouch for use with breaching tools or storage for a hydration pack
- Expansion zipper to increase load capacity
- External MOLLE on sides and bottom of bag / Internal MOLLE rows on inside back of bag
- Three clear tear away pouches on the inside main compartment
- Two internal mesh pouches and one larger internal pouch on inner lid
- Internal tape roll strap

Inside Back Compartment for stowing extra equipment



Every MINI MEDIC KIT includes these core items:



SPECIFICATIONS

H 14 in. x W 11 in. x D 5 in.
(Depth expandable to 7.25 In.)
Weight: 12.5 lb

MINI MEDIC KIT

COLOR	LEVEL	ITEM #
● Black	Basic	80-0979
● Coyote	Basic	80-0980
● OD Green	Basic	80-0981
● Ranger Green	Basic	80-1003
● MultiCam	Basic	80-1231
● Black	Combat Gauze	80-0985
● Coyote	Combat Gauze	80-0986
● OD Green	Combat Gauze	80-0987
● Ranger Green	Combat Gauze	80-1005
● MultiCam	Combat Gauze	80-1232
● Black	Advanced	80-0982
● Coyote	Advanced	80-0983
● OD Green	Advanced	80-0984
● Ranger Green	Advanced	80-1004
● MultiCam	Advanced	80-1233

Bag only option also available at NARescue.com

PLUS+ BASIC
also includes:



COMBAT GAUZE
also includes:



ADVANCED also includes:



Patrol **VEHICLE TRAUMA KIT**

- Provide personnel with an easily accessible, quick to deploy medical kit from your Patrol Vehicle
- Modular storage system comes organized with an array of medical products ranging from basic first aid to more intermediate trauma care
- One side covers all your agency needs to administer basic first aid and community assistance
- Second side offers a first time combination of Basic First Aid and Basic Trauma medical items
- Exterior includes MOLLE-style system to accommodate extra equipment as needed



PATROL VEHICLE TRAUMA KIT:

 QTY 1 C-A-T® Tourniquet	 QTY 1 ETD™ 6 in. Responder	 QTY 1 ETD™ Abdominal	BASIC QTY 2 Gauze Compressed	-OR- COMBAT GAUZE QTY 1 Gauze Compressed	 QTY 1 Combat Gauze® LE
 QTY 1 CPR Mask	 QTY 1 HyFin® Vent, Twin Pk	 QTY 1 Survival Blanket	 QTY 1 BurnTec 4"x4" Dressing	 QTY 4 Gauze Pads 4"x4" (2pk)	 QTY 25 Flexible Fabric Bandages
 QTY 2 Triangular Bandage	 QTY 2 Elastic Bandage, 4 in.	 QTY 1 SAM® Splint II	 QTY 1 Trauma Shears, Lg	 QTY 2 Surgical Tape, 2" Mini	 QTY 1 Eye Wash, 4 oz
 QTY 2 PES Eye Shield	 PAIR 3 Nitrile Gloves, Lg	 QTY 1 Pen Light	 QTY 1 Forceps, Splinter, 4.5 in.	 QTY 1 Permanent Marker, Large	

PATROL VEHICLE TRAUMA

LEVEL	ITEM #
Basic	80-0542
Combat Gauze LE	80-0859

SPECIFICATIONS

H 7 in. x W 12 in. x D 5 in.
Weight: 6 lb



Trauma Kits

RANGE TRAUMA AID KIT Hard Case

- Water-tight case maintains the integrity of your medical products when you need them most
- Kit contents arranged in a modular system featuring rugged lid storage organizer and removable nylon pouches
- Labeled internal pouches for ease of use



RANGE TRAUMA AID

ITEM #

85-0889

SPECIFICATIONS

H 9.5 in. x W 13.5 in. x D 6 in.

Weight: 8 lb 12 oz

RANGE TRAUMA AID KIT HARD CASE:



School Resource Officer Kits

SRO CRISIS RESPONSE KIT

- Equip School Resource Officers with the essential life-saving tools to begin treating injured victims until medical personnel arrive
- Empowers an SRO (School Resource Officer) to administer self-aid in the event he or she is injured or provide immediate care to injured casualties during a school incident
- Multiple carry options

SRO CRISIS RESPONSE KIT:



Bag only option also available at NARescue.com



SPECIFICATIONS

H 7.75 in. x W 13 in. x D 9 in.

Weight: 5 lb 16 oz

SRO CRISIS RESPONSE KIT

ITEM #

85-0417

I-ROK *Individual Rugged Outdoor Kit*



- Contains essential first-aid supplies to treat basic injuries as well as control major bleeding
- Storage case is crushproof, dustproof and waterproof up to a depth of 2 meters
- Flexible lid organizer helps keep smaller medical supplies neatly stored
- Built-in pressure valve keeps water and dust out while balancing air pressure
- Available in 4 hardcase colors and 2 sizes with the option of Basic or with Bleeding Control Dressing (hemostatic bandage) versions

Every I-ROK contains the following items:



QTY 1 C-A-T® Tourniquet

BASIC QTY 1 Responder Compressed Gauze

-OR- BASIC with Bleeding Control Dressing QTY 1 Bleeding Control Dressing



I-ROK

COLOR	LEVEL	ITEM #
Black	Basic	82-0031
Coyote	Basic	82-0032
OD Green	Basic	82-0033
Orange	Basic	82-0034
Black	with BCD	82-0035
Coyote	with BCD	82-0036
OD Green	with BCD	82-0037
Orange	with BCD	82-0038

SPECIFICATIONS

H 7 in. x W 4.5 in. x D 2.5 in.
Weight: 2.48 lb



QTY 1 ETD™ 4 in. Mini

QTY 1 HyFin® Vent Compact, Twin

QTY 5 Butterfly Bandage

QTY 5 1x3 in. Bandage

QTY 3 Burn Gel



QTY 1 MicroShield CPR Face Mask

QTY 3 Antibiotic Ointment

QTY 2 Sting & Bite Packet

QTY 1 Survival Blanket

QTY 3 Hand Sanitizer

PAIR 1 Responder Nitrile Gloves, Lg

QTY 1 Mini Duct Tape, 2 in.

QTY 1 Permanent Marker, Large

QTY 1 Just in Time Instructions

I-ROK XL

SPECIFICATIONS

H 8.9 in. x W 5.4 in. x D 3 in.
Weight: 3 lb

I-ROK XL includes these extra items:



QTY 5 Butterfly Bandage

QTY 5 1x3 in. Bandage

QTY 2 Burn Gel



QTY 1 BurnTec 4"x4" Dressing

QTY 2 Antibiotic Ointment

QTY 1 Sting & Bite Packet

QTY 2 Hand Sanitizer

QTY 2 Gauze Pads 4"x4" (2pk)

QTY 1 Triangular Bandage

PAIR 1 Responder Nitrile Gloves, Lg

QTY 1 Forceps, Splinter, 4.5 in.

QTY 1 Trauma Shears, Lg

I-ROK XL

COLOR	LEVEL	ITEM #
Black	Basic	82-0049
Coyote	Basic	82-0050
OD Green	Basic	82-0051
Orange	Basic	82-0052
Black	with BCD	82-0053
Coyote	with BCD	82-0054
OD Green	with BCD	82-0055
Orange	with BCD	82-0056

Rescue Tools & Litters



Leatherman® RAPTORM™ MULTI-TOOL

High-performance multi-tool equipped with a variety of mission critical tools essential to emergency medical practitioners in the field

1. 420HC Stainless Steel Folding Medical Shears
2. Strap Cutter
3. Oxygen Tank Wrench
4. Ring Cutter
5. Carbide Glass Breaker
6. 5 cm Ruler
7. Holster with Clip

RAPTOR MULTI-TOOL

HANDLE	ITEM #	NSN #
● Black	70-0014	5110-01-626-8766
● Blue	70-0264	----
● Coyote	70-0265	----
● Red	70-0266	----
● Green	70-0326	----
● Org/Blk	70-0327	----

SPECS


Shear Length: 1.9 in.
Closed Length: 5 in.
Weight: 5.8 oz

Leatherman® RAPTORM™ RESPONSE SHEARS

1. Stainless Steel folding shears
2. Integrated ring cutter
3. 5 cm ruler
4. Oxygen tank wrench
5. Lanyard connection in handle
6. Removable pocket clip

RAPTOR RESPONSE SHEARS

HANDLE	ITEM #
● Cement	70-0328
● Crimson	70-0329
● Navy	70-0330



SPECS

Length: Closed 4.3 in.
Open 7.25 in.
Shear Length: 1.9 in.
Width: 1.7 in | 4.3 cm
Overall Thickness: .85 in.
Weight: 5.5 oz

Benchmade® HOOK SAFETY CUTTER

- Saves time – more than 1.5 minutes faster than traditional trauma shears when cutting 2 articles of clothing.
- Extensively tested and widely used by the US military for medical and egress purposes
- Maintains sharpness for hundreds of cuts (up to 500+ cuts through seatbelt test material)
- 440C cutting surface and opening with handle to provide a more compact carry
- Includes a MOLLE-compatible carrying sheath



HOOK SAFETY CUTTER

ITEM #	STYLE
ZZ-0652	7 Hook

SPECS

Length: 4.3 in.
Weight: 1.5 oz

NAR QUIKLITTER LITE

- Perfect balance of size, weight, cost, and capability
- Features ten large carry handles with a load capacity tested to 750 lb
- Made of lightweight, heavy-duty nylon
- Vacuum-packed to conserve space and fit in a cargo pocket
- Latex-free & water-resistant

QUIKLITTER

SIZE	ITEM #
LITE*	60-0071

*Nylon Construction

SPECIFICATIONS

Payload 750 lb.
Litter Bed: L 74 in. x W 27 in.
Weight: 12 oz



UNO EMERGENCY SPLINT

- Modular Design
- Stepless Adjustment
- Small and Lightweight
- Compatible with TCCC Equipment
- Pack of 3 includes additional ankle strap for stabilizing the ankle and enables the device to provide traction



SPECIFICATIONS

Packaged - L 21.50 in. x W 10.55 in. x D 1.50 in.
Unpackaged - L 15.75 in. x W 2.17 in. x D 2.36 in.
Weight: 0.91 lbs

UNO EMERGENCY SPLINT

ITEM #	OPTION
50-0045	Single
50-0046	Pack of 3



NORTH AMERICAN RESCUE®

Education and Training


North American Rescue Education and Training (NAR Training)

has roots that run deep and wide in the world of pre-hospital medicine and technical rescue. Our team was built from the ground up on the vision and experience of seasoned military medics seeking to fill the training gaps commonly encountered in the operational environment.

Since then, we have continued to broaden our vision to include other invaluable first responder perspectives such as pre-hospital healthcare providers in civilian, law enforcement, emergency medical services and community public safety initiatives.

CORE Competencies

- Continuing Medical Education (CMEs)
- Prolonged Casualty Care
- Austere Medicine
- Technical Rescue
- Emergency Medical Services
- Critical Care Flight Medicine



Training Locations:

- NAR Training
- Training Partners & Affiliates

ACCREDITATIONS:



American Heart Association®


Now SAME DAY Cards Issued At End Of Course!



Explore training opportunities online

NAEducation.com ➤ **Phone: 888.689.6277** ➤ **Email: NARTraining@NARescue.com**

CAT[®] COMBAT APPLICATION TOURNIQUET



**THE FASTEST,
SAFEST, MOST
EFFECTIVE
PREHOSPITAL FIELD
TOURNIQUET**

TOURNIQUETS | FIRST AID KITS | TRAUMA KITS | RESCUE TOOLS | LITTERS



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

35 Tedwall Court | Greer, SC 29650-4791