

PUBLIC ACCESS  
**BLEEDING CONTROL**

INDIVIDUAL KITS | PACKS | STATIONS



PUBLIC ACCESS BLEEDING CONTROL KITS | BLEEDING CONTROL TRAINERS

# Every Second Counts...



COMMUNITY  
PREPAREDNESS



## 3 to 5 MINUTES

### *The Difference between Life and Death*

- ▶ A victim who is bleeding can quickly die BEFORE professional emergency responders can even arrive on the scene
- ▶ Serious bleeding accounts for 35% of prehospital deaths from traumatic injuries\*
- ▶ Life-threatening bleeding warrants immediate intervention
- ▶ Basic bleeding control provided by bystanders can often be the difference between life and death...



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

Trauma is a national health problem and after a traumatic injury, serious bleeding is responsible for over 35% of prehospital deaths\*. Examples of trauma include traffic & industrial accidents, natural disasters, falls, accidental injuries and intentional acts of violence such as shootings, stabbings and bombings.

These types of injuries often result in uncontrolled bleeding. Without immediate intervention these injuries can result in death within minutes – even before emergency medical services or other first responders can arrive.

### **North American Rescue**

is leading the way to equip bystanders and initial first responders with the essential medical gear needed to take immediate action in life threatening emergencies. These lifesaving modular aid components are designed for individuals or groups in the workplace, public buildings, and populated public areas to administer immediate bleeding control at the point of wounding.

**Death from Uncontrolled BLEEDING can occur in as little as 3 to 5 Minutes**

\*Chambers JA, Seastedt K, Krell R, Caterson E, Levy M, Turner N. "Stop the Bleed": A U.S. Military Installation's Model for Implementation of a Rapid Hemorrhage Control Program. Mil Med. 2019 Mar 01;184(3-4):67-71. [PubMed] [Reference list] <https://academic.oup.com/milmed/article/184/3-4/67/5062295>

# Individual Kits | Public Access Bleeding Control



PVC  
**BLEEDING CONTROL PATCH**

- Soft flexible PVC rubber patch for labeling your lifesaving medical gear
- Hook-&-loop attachment

**SPECS**

L 4.375 in. x H 1.75 in.  
x D 0.125 in.  
Weight: 0.275 oz

**PATCH**

**ITEM #**  
ZZ-0796

## Public Access **INDIVIDUAL BLEEDING CONTROL KITS**

### Nylon Bag Option:

- Easy to restock in cases of use
- Convenient to carry and easy to open in stressful situations
- Opens flat for ease of use and immediate identification of enclosed devices
- Easily recognizable red nylon bag with attached PVC Bleeding Control ID Patch

### Vacuum Sealed Option:

- Cost-effective alternative to Nylon
- Compact and rugged Vacuum-Sealed pouch
- Large Identifying Card Vacuum-Sealed inside the pouch for ease of identification



**SPECIFICATIONS**

Nylon Bag: H 7 in. x W 5 in. x D 2.5 in.  
Vacuum Sealed: H 7 in. x W 5 in. x D 2.75 in.

**PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KITS**

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	15 oz	80-0452
Intermediate	Red Nylon	16.5 oz	80-0453
Advanced BCD	Red Nylon	17.4 oz	80-0523
Advanced	Red Nylon	17.4 oz	80-0454
Basic	Vacuum Sealed	11.7 oz	80-0465
Intermediate	Vacuum Sealed	13.25 oz	80-0466
Advanced BCD	Vacuum Sealed	14 oz	80-0522
Advanced	Vacuum Sealed	14 oz	80-0467

Every **INDIVIDUAL BLEEDING CONTROL KIT** contains:



**PLUS** these Individual Components based upon YOUR configuration choice:



## Trauma & First Aid **ESSENTIAL BLEEDING CONTROL KIT**

Every **ESSENTIAL BLEEDING CONTROL KIT** contains:



**PLUS+** one of these:



- Compact and Convenient Kit
- Highly mobile and easily deployed
- Hypothermia prevention included
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the safest and most effective Tourniquet on the market



**SPECIFICATIONS**

L 12 in. x H 7.5 in.  
Weight: 0.68 lbs

**TRAUMA & FIRST AID ESSENTIAL BLEEDING CONTROL KIT**

OPTION	ITEM #
Compressed Gauze	80-1474
Bleeding Control Dressing	80-1475
CELOX™ First Aid Gauze	80-1476

Public Access

# BLEEDING CONTROL MULTI-PACKS



## 5-PACKS *Vacuum*

- Five (5) Vacuum-sealed Individual Bleeding Control Kits inside an easy-to-open nylon tote
- Includes One (1) Responder QuikLitter™ for moving patients
- Available in Basic, Intermediate, Advanced BCD & Advanced



### SPECIFICATIONS

#### 5-Pack:

H 7.5 in. x W 13 in. x D 8 in.

Each Individual

Bleeding Control Kit:

H 7 in. x W 5 in. x D 2.5 in

### PUBLIC ACCESS BLEEDING CONTROL 5-PACKS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Vacuum Sealed	5.7 lb	80-0573
Intermediate	Vacuum Sealed	5.75 lb	80-0574
Advanced BCD	Vacuum Sealed	6.5 lb	80-0576
Advanced	Vacuum Sealed	6.5 lb	80-0575

## 8-PACKS *Nylon | Vacuum Sealed*

- Eight (8) Nylon or Vacuum-sealed Individual Bleeding Control Kits inside an easy-to-open nylon carrying case
- Includes One (1) Responder QuikLitter™ for moving patients
- Available in Basic, Intermediate, Advanced & Advanced BCD

### SPECIFICATIONS

#### 8-Pack:

H 14 in. x W 11 in. x D 8 in.

Each Individual Bleeding Control Kit:

H 7 in. x W 5 in. x D 2.5 in.

### PUBLIC ACCESS BLEEDING CONTROL 8-PACKS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	8.5 lb	80-0460
Intermediate	Red Nylon	9.3 lb	80-0461
Advanced BCD	Red Nylon	9.7 lb	80-0525
Advanced	Red Nylon	9.7 lb	80-0462
Basic	Vacuum Sealed	6.8 lb	80-0468
Intermediate	Vacuum Sealed	7.6 lb	80-0469
Advanced BCD	Vacuum Sealed	8 lb	80-0524
Advanced	Vacuum Sealed	8 lb	80-0470

## TWIN PACKS

- Contents of 2 Individual Kits
- Compact Nylon Bag designed to fit next to an AED in the cabinet



### SPECIFICATIONS


Nylon Bag:

H 10 in. x W 9 in. x D 3 in.

### TWIN PACKS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	1.66 lb	80-0951
Intermediate	Red Nylon	1.8 lb	80-0952
Advanced BCD	Red Nylon	1.8 lb	80-0953

Every **BLEEDING CONTROL TWIN PACK** contains:



**PLUS** these Individual Components based upon YOUR configuration choice:



# Clear Polycarbonate WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in red nylon bags or Vacuum-Sealed pouches.
- The 5 & 8 Pack Stations also include an easy-to-open nylon carrier and one (1) Responder QuikLitter™ for moving injured casualties.

## SPECIFICATIONS

**5 & 8 Pack Stations:** H 18 in. x W 15.5 in. X D 10 in.  
Empty Case Wt: 15 lb (includes rack)  
Large Case Item# ZZ-0822

**Individual Stations:** H 10 in. x W 10 in. X D 5.25 in.  
Empty Case Wt: 2.39 lb  
Small Case Item# ZZ-1000

## 8-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	21.8 lb	80-0471
Intermediate	8	Vacuum	22.6 lb	80-0472
Advanced BCD	8	Vacuum	23 lb	80-0526
Advanced	8	Vacuum	23 lb	80-0473
Basic	8	Nylon	23.5 lb	80-0474
Intermediate	8	Nylon	24.3 lb	80-0475
Advanced BCD	8	Nylon	27.4 lb	80-0527
Advanced	8	Nylon	27.4 lb	80-0476

## 5-PACK CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	5	Vacuum	20.7 lb	80-0581
Intermediate	5	Vacuum	20.7 lb	80-0582
Advanced BCD	5	Vacuum	21.5 lb	80-0584
Advanced	5	Vacuum	21.5 lb	80-0583

## IND. CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	1	Vacuum	3.3 lb	80-0938
Intermediate	1	Vacuum	3.42 lb	80-0939
Advanced BCD	1	Vacuum	3.48 lb	80-0940
Basic	1	Nylon	3.3 lb	80-0935
Intermediate	1	Nylon	3.6 lb	80-0936
Advanced BCD	1	Nylon	3.7 lb	80-0937

**Alarm (#ZZ-1075) is also available for the Clear Polycarbonate Wall Cabinet**

Cases are water resistant. To best protect your critical lifesaving products, cases are best mounted under cover and out of the direct elements.



8 Pack Station



Individual Kit Station



5 Pack Station



Small Individual Clear Case: ZZ-1000



Large Multi-Pack Clear Case: ZZ-0822



Designed to be mounted on a wall at a designated location or to be co-located with Automatic External Defibrillators (AEDs) to provide publically accessible medical supplies to control major bleeding

- Durable, clear polycarbonate cabinet with rounded corners provides both physical and environmental protection for medical equipment
- Large case features wire shelf to support the contents
- External eyelet allows a secondary safety seal to easily identify if the station has been opened or tampered with
- External mounting tabs and hardware enable easy installation on any wall
- Clear design allows visibility at any angle

**CLEAR CASE features**



5 Pack Station

**BLEEDING CONTROL KIT**



## Semi-Recessed Metal WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in either red nylon bags or Vacuum-Sealed pouches stored within an easy to-open nylon carrier and one (1) Responder QuikLitter™ for moving injured patients
- Outside Dimensions: H 17.5 in. x W 17.5 in. x D 1 in.  
Inside Dimensions: H 14.13 in. x W 14.13 in. x D 6.88 in.  
Empty Case Wt: 10 lb (Empty Case Item# ZZ-0921)



8 Pack Station

**METAL CASE**  
features

- Mounts permanently inside the wall with a low profile 1 in. trim
- Aluminum cabinet with a white finish and clear plexiglass window
- Stainless steel continuous hinge & chrome plated handle designed to open with one quick pull

**8-PACK SEMI-RECESSED METAL WALL STATIONS**

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	16.8 lb	80-0861
Intermediate	8	Vacuum	17.6 lb	80-0862
Advanced BCD	8	Vacuum	18 lb	80-0863
Advanced	8	Vacuum	18 lb	80-0864
Basic	8	Nylon	18.5 lb	80-0865
Intermediate	8	Nylon	19.3 lb	80-0866
Advanced BCD	8	Nylon	15.7 lb	80-0867
Advanced	8	Nylon	15.7 lb	80-0868

**5-PACK SEMI-RECESSED METAL WALL STATIONS**

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	5	Vacuum	15.7 lb	80-0577
Intermediate	5	Vacuum	15.7 lb	80-0578
Advanced BCD	5	Vacuum	16.5 lb	80-0580
Advanced	5	Vacuum	16.5 lb	80-0579



**Semi-Recessed Metal Case: ZZ-0921**

*All Multi-Pack Wall Stations feature: Highly visible "Bleeding Control Kit" letters for fast ID Instructions for Use posted on the interior back panel*

*Cases are water resistant. To best protect your critical lifesaving products, cases are best mounted under cover and out of the direct elements.*



**8 Pack Nylon Station**



**Low Profile Metal Cabinet: ZZ-1199**

*Designed to be mounted on a wall with a slim 4 in. clearance and featuring easy pull-to-open design for rapid public access to major bleeding control supplies*

## Low Profile Metal WALL STATIONS

- Low profile metal cabinet mounts on the wall with a slim 4 in. clearance
- Designed to safely store 8 Individual Bleeding Control Kits in public areas
- Kits are available in red nylon bags or Vacuum-Sealed pouches and stored in Two (2) easily removable Mesh Panels each equipped with 4 pockets and nylon handle
- Features highly visible "Bleeding Control Kit" lettering for fast identification and convenient Instructions for Use on the interior door panel

### SPECIFICATIONS

H 15.5 in. x W 22.5 in. X D 4 in.

Empty Cabinet Weight: 9.365 lb

**Low Pro Cabinet with ID Stickers:**  
Item# ZZ-1199



**8 Pack Vacuum Sealed Station**

### 8-PACK LOW PROFILE METAL WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	17.49 lb	80-1083
Intermediate	8	Vacuum	18.19 lb	80-1084
Advanced BCD	8	Vacuum	18.69 lb	80-1085
Basic	8	Nylon	17.49 lb	80-1080
Intermediate	8	Nylon	18.19 lb	80-1081
Advanced BCD	8	Nylon	18.69 lb	80-1082

**Alarm (#ZZ-1075) is also available for the Low Profile Metal Cabinet**

## Low-Pro CABINET features

- Low profile metal cabinet with 4 in. wall clearance
- Provides convenient slimline storage for medical equipment
- White powder painted corrosive resistant galvanized steel cabinet with stainless steel continuous hinge
- Magnetic closure provides quick easy access with one pull to open - No twist and pull

Cases are water resistant. To best protect your critical lifesaving products, cases are best mounted under cover and out of the direct elements.



- Small, compact and easy to use
- Kits contain the same Tourniquet and Hemostatic Dressing technology recommended by the U.S. Military & the American College of Surgeons Committee on Trauma to control bleeding
- Available in Vacuum-sealed Pouch or Red Nylon Bag

## STB KITS

Nylon or Vacuum-Sealed



### VAC STB KIT

ITEM #

85-1574

### SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.

Weight: 8 oz

#### VACUUM-SEALED STB KIT contains:



QTY  
**1** C-A-T®  
Tourniquet

QTY  
**1** ETD™ 4 in.  
Responder, Flat

QTY  
**1** Bleeding  
Control Dressing

PAIR  
**1** Responder  
Nitrile Gloves, Lg

QTY  
**1** Just in Time  
Instructions

QTY  
**1** Permanent  
Marker, Small



### SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.

Weight: 15 oz

### NYLON STB KIT

ITEM #

85-2995

#### NYLON STB KIT contains:



QTY  
**1** C-A-T®  
Tourniquet

QTY  
**1** Responder  
6 in. ETD

QTY  
**1** Wound  
Packing Gauze

QTY  
**1** Bleeding  
Control Dressing

PAIR  
**2** Responder  
Nitrile Gloves, Lg

QTY  
**1** Just in Time  
Instructions

QTY  
**1** Permanent  
Marker, Small

## Cabinet Security ALARMS

- Alarm compatible with Clear Polycarbonate and Low Profile Metal Cabinets
- Loud 120 decibel alarm sounds when magnet separates from base more than 0.5 in.
- Easy and convenient to install
- Requires no wiring
- Mounting hardware and instructions included
- Designed for indoor use only



### SPECIFICATIONS

W 1.63 in. x H 2.41 in. x D 0.72 in.

Weight: 6.4 oz each

### CABINET SECURITY ALARMS

ITEM#	UNIT
ZZ-1075	Pack of 2

## Accessories

### 3-WAY SIGN Bleeding Control



- Highly visible for locating lifesaving bleeding control kits immediately
- 3 way "V" shape projects from wall for easy identification from multiple directions

### SPECIFICATIONS

H 7 in. x W 16.5 in. x D 6.75 in.

Weight: 0.35 oz

### 3-WAY SIGN

ITEM #

ZZ-1117

# Workplace First Aid Kits

## Class A **TRAUMA & FIRST AID KIT**

- Designed as a portable First Aid Kit for the workplace **exceeding ANSI/ISEA Z308.1-2021 Standard for Class A kits**
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the C-A-T Tourniquet - the safest and most effective Tourniquet on the market!
- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side



### TRAUMA & FIRST AID KIT - CLASS A

CASE	DRESSING	ITEM #
Soft Case	Compressed Gauze	80-0947
Soft Case	Bleeding Control Dressing	80-1045
Soft Case	CELOX™ First Aid Gauze	80-1477
Hard Case	Compressed Gauze	80-1033
Hard Case	Bleeding Control Dressing	80-1046
Hard Case	CELOX™ First Aid Gauze	80-1481

### SPECIFICATIONS

**Soft Case:**  
L 14.75 in. x H 8 in. x D 4 in.  
Weight: 2.94 lb

**Hard Case:**  
L 13 in. x H 9.5 in. x D 5 in.  
Weight: 6 lb

Type IV ANSI Container only option (#80-1113) also available at [NARescue.com](http://NARescue.com)

### TRAUMA & FIRST AID KIT CLASS A:

 QTY 1 C-A-T® Tourniquet	 QTY 1 HyFin® Vent Compact Twin	 QTY 1 Responder Compressed Gauze	 QTY 1 ETD™ 4 in. Responder, Flat	 QTY 1 4" Elastic Bandage	 QTY 1 2" Elastic Bandage	 QTY 1 Triangular Bandage
 QTY 1 SAM® Splint II	 QTY 2 4" x 4" Gauze Pad (2pk)	 QTY 2 Gauze Pad, 5" x 9"	 QTY 16 1" x 3" Bandage	 QTY 1 BurnTec 4"x4" Dressing	 QTY 1 CPR Microshield Mask	 QTY 10 Antiseptic Towelette
 QTY 10 Triple Antibiotic Ointment	 QTY 1 Surgical Tape, 1 in.	 QTY 1 6"x9" Instant Cold Pack	 QTY 1 Eye Wash Solution - 1 oz	 QTY 2 Sterile Oval Eye Pad	 PAIR 2 Nitrile Gloves, Lg	 QTY 1 Responder Shears, Lg
 QTY 10 Hand Sanitizer	 QTY 10 Safetec Burn Gel	 QTY 1 SURVIVAL BLANKET	 QTY 1 Permanent Marker, Large	 QTY 1 First Aid Pocket Guide	 QTY 1 INSTRUCTIONS FOR USE PUBLIC BLEEDING CONTROL	 QTY 1 TRAUMA & FIRST AID KIT - CLASS A KIT CONTENTS



**PLUS+**  
one of these:

 w/ Comp. Gauze QTY 1 Responder Compressed Gauze	- OR -	 w/ BCD QTY 1 Bleeding Control Dressing	- OR -	 w/ CELOX™ QTY 1 5ft CELOX Hemostatic Dressing
--	--------	---	--------	--

### TRAUMA & FIRST AID KIT RESUPPLY - CLASS A

DRESSING	ITEM #
Compressed Gauze	80-1153
Bleeding Control Dressing	80-1155
CELOX™ First Aid Gauze	80-1479

Container not included with ReSupply Kits

# Workplace First Aid Kits

## Class B **TRAUMA & FIRST AID KIT**

- Designed as a portable First Aid Kit for the workplace **exceeding ANSI/ISEA Z308.1-2021 Standard for Class B kits**
- Same supplies as the Class A Kit but with larger quantities to accommodate larger workplaces
- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side

### SPECIFICATIONS

#### Soft Case:

L 14.75 in. x H 8 in. x D 7.5 in.  
Weight: 4.87 lb

#### Hard Case:

L 13.5 in. x H 11 in. x D 6 in.  
Weight: 9 lb

### TRAUMA & FIRST AID KIT - CLASS B

CASE	DRESSING	ITEM #
Soft Case	Compressed Gauze	80-0948
Soft Case	Bleeding Control Dressing	80-1047
Soft Case	CELOX™ First Aid Gauze	80-1478
Hard Case	Compressed Gauze	80-1034
Hard Case	Bleeding Control Dressing	80-1048
Hard Case	CELOX™ First Aid Gauze	80-1482

Type IV ANSI Container only option (#80-1114) also available at [NARescue.com](http://NARescue.com)

### TRAUMA & FIRST AID KIT CLASS B:

 QTY 1 C-A-T™ Tourniquet	 QTY 1 HyFin® Vent Compact Twin	 QTY 1 Responder Compressed Gauze	 QTY 1 ETD™ 4 in. Responder, Flat	 QTY 1 4" Elastic Bandage	 QTY 2 2" Elastic Bandage	 QTY 2 Triangular Bandage	
 QTY 1 SAM® Splint II	 QTY 4 4" x 4" Gauze Pad (2pk)	 QTY 4 Gauze Pad, 5" x 9"	 QTY 50 1" x 3" Bandage	 QTY 2 BurnTec 4" x 4" Dressing	 QTY 1 CPR Microshield Mask	 QTY 50 Antiseptic Towelette	
 QTY 25 Triple Antibiotic Ointment	 QTY 2 Surgical Tape, 1 in.	 QTY 2 6" x 9" Instant Cold Pack	 QTY 1 Eye Wash Solution - 1 oz	 QTY 2 Sterile Oval Eye Pad	 PAIR 4 Responder Nitrile Gloves, Lg	 QTY 1 Responder Shears, Lg	
 QTY 20 Hand Sanitizer	 QTY 25 Safetec Burn Gel	 QTY 1 Survival Blanket	 QTY 1 Permanent Marker, Large	 QTY 1 First Aid Pocket Guide	 QTY 1 INSTRUCTIONS FOR USE PUBLIC BLEEDING CONTROL	 QTY 1 TRAUMA & FIRST AID KIT - CLASS A Kit Contents Inventory Card	 QTY 1 Kit Contents Inventory Card

### TRAUMA & FIRST AID KIT RESUPPLY - CLASS B

DRESSING	ITEM #
Compressed Gauze	80-1154
Bleeding Control Dressing	80-1156
CELOX™ First Aid Gauze	80-1480

Container not included with ReSupply Kits

### PLUS+ one of these:

 w/ Comp. Gauze QTY 1 Responder Compressed Gauze	- OR -	 w/ BCD QTY 1 Bleeding Control Dressing	- OR -	 w/ CELOX™ QTY 1 5ft CELOX Hemostatic Dressing
--	--------	--	--------	---

# Public Access INDIVIDUAL BLEEDING CONTROL TRAINER KITS

- The appropriate equipment to teach individuals the necessary life-saving skills of bleeding control
- Corresponds with NAR's single Public Access Bleeding Control Kits in Basic, Intermediate, or Advanced
- Packaged in a compact, blue nylon bag with PVC Bleeding Control Patch

Every **INDIVIDUAL TRAINER KIT** contains:



Bag only option available at [NARescue.com](http://NARescue.com)

### SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.

PLUS these Individual Components based upon YOUR configuration choice:



### BLEEDING CONTROL INDIVIDUAL TRAINER KITS

LEVEL	WGT	ITEM #
Basic	15 oz	80-0901
Intermediate	16.5 oz	80-0902
Advanced	16.5 oz	80-0903



## STB TRAINER

- Designed to provide the essential training equipment to teach how to control bleeding
- Contains training items corresponding with the necessary tools recommended by the American College of Surgeons Committee on Trauma to equip bystanders to help save lives in the instance of life threatening bleeding

**STB TRAINER** contains:



### SPECIFICATIONS

H 7 in. x W 5 in. x D 2.5 in.

Weight: 8 oz

### STB TRAINER

ITEM #
85-1880

# Multi-User Skills Training Kits | Public Access Bleeding Control

## BASIC BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



## INTERMEDIATE BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



## ADVANCED BLEEDING CONTROL TRAINING SKILLS KIT CONTENTS:



# BLEEDING CONTROL TRAINING SKILLS KITS



- Provides the appropriate equipment for up to 20 students to learn the necessary lifesaving skills of bleeding control
- Packaged in a self-contained, wheeled, rigid case with internal dividers for easy organization and inventory
- Skill set training addresses extremity bleeding, wound packing, pressure bandages, chest injuries and patient transport
- Corresponds with NAR's Public Access Bleeding Control Kit Levels: Basic, Intermediate & Advanced
- Available with or without Wound Mgt. Simulator

### SPECIFICATIONS

Hard Storage Case:  
W 32 in.x H 21 in. x D 12.5 in.

### BLEEDING CONTROL TRAINING SKILLS KITS

LEVEL	WGT	ITEM #
Basic	40 lb	80-0881
Intermediate	42 lb	80-0879
Advanced	42 lb	80-0877
Basic   Wound Mgt. Simulator	44 lb	80-0880
Intermediate   Wound Mgt. Simulator	48 lb	80-0878
Advanced   Wound Mgt. Simulator	48 lb	80-0876



## Bleeding Control Training Aids

# TRUECLOT® Hemorrhage Training Kit™

- Introduce students to realistic hemorrhage control for TCCC, TECC, TEMS, PHTLS, Active Shooter Training and general wound packing techniques for advanced hemostatic gauze training
- Highly realistic Wound Packing Trainer features a large caliber gunshot wound with hidden bleed points and a simulated fracture/defect hard plastic bone insert to aid in teaching deep wound packing techniques
- Includes TrueClot® Clotting Blood Simulant and Simulated Hemostatic Gauze which demonstrates the reaction of Hemostatic Gauze in a wound
- Durable design features a sturdy silicone construction with no moving parts, reinforced tubing mount, and a more realistic wound pattern and outer anatomy, ideal for the repeated use of a training environment

### TRUECLOT HEMORRHAGE TRAINING KIT

#### ITEM #

93-0020

#### SPECIFICATIONS

Task Trainer: L 7.25 in. x H 2.5 in.

Liquid Capacity: 1 liter

Storage Case: L 13.5 in. x W 10 in. x D 5.5 in.

Weight: 5 lb



### TRUECLOT TRAINER:



# TRUECLOT® Tourniquet Simulator

- Wearable device that allows full application of a tourniquet for training purposes
- Ideal for training tourniquet placement techniques featuring a durable silicone construction with no moving parts, reinforced tubing mount, and a more realistic wound pattern and outer anatomy

### TRUECLOT TOURNIQUET SIMULATOR



#### SPECIFICATIONS

H 12 in. x W 16 in. x D 7 in.

Weight: 5 lb

### TRUECLOT TQ SIMULATOR

#### ITEM #

93-0931



## TRUECLOT® Bleeding Control Instructor Simulator Kit

with 1 Gun Shot Wound  
& 1 Knife/Laceration

- Everything needed to conduct a wet or dry wound packing training class
- Designed to simulate wound packing and demonstrate the use of advanced hemostatic gauzes with LUNA's TrueClot® Blood Simulant
- Task Trainer for teaching realistic hemorrhage control in wounds including a large caliber gunshot & knife/laceration
- Training capabilities include general wound packing techniques & advanced hemostatic gauze training

### TRUECLOT BLEEDING CONTROL INSTRUCTOR SIMULATOR

ITEM #

93-0932



### SPECIFICATIONS

H 17 in. x W 20 in. x D 7.5 in.  
Weight: 14 lb

### TRUECLOT GSW & Laceration SIMULATOR:



QTY  
**1** GSW  
Wound Model



QTY  
**1** Laceration  
Wound Model



QTY  
**4** TrueClot®  
Blood Simulant



QTY  
**10** TrueClot®  
Hemostatic, 12 ft




QTY  
**2** Wound  
Pans  
for Blood  
Containment



QTY  
**5** TrueClot®  
Flow Valves



QTY  
**2** MedMat  
Sterile Procedural Mat



QTY  
**2** Bags,  
Disposable



QTY  
**2** TrueClot®  
Mixing Bottles



QTY  
**2** TrueClot®  
Bleeding Bottles



QTY  
**8** Absorbent  
Underpad Chux



QTY  
**1** Storage &  
Transport Case



QTY  
**20** Nitrile  
Gloves (Medium)



QTY  
**20** Nitrile  
Gloves (Large)



QTY  
**20** Nitrile  
Gloves (X-Large)

## Bleeding Control Training Aids



### COMBAT GAUZE<sup>®</sup> Moulage Trainer

- Thoroughly instruct students on the proper application of Z-Fold Combat Gauze<sup>®</sup>
- Completely inert (no hemostatic properties) - training purposes only
- Re-sealable packaging for multiple uses

#### SPECIFICATIONS


Packaged: H 6 in. x W 4.5 in. x D 1.25 in.  
Gauze: L 4 yd x W 3 in.  
Weight: 0.75 oz

#### COMBAT GAUZE TRAINER

##### ITEM #

30-0063

### C-A-T<sup>®</sup> TRAINER Combat Application Tourniquet



- Fully functioning C-A-T<sup>®</sup> tourniquet colored blue to identify it as a training device
- Enables instructors and trainees hands-on training & tourniquet application practice
- Providing a distinct visual means to distinguish training devices from personal/operational life-saving devices

#### C-A-T<sup>®</sup> TRAINER

##### ITEM #

30-0033

#### SPECIFICATIONS

H 6.5 in. x W 2.4 in. x D 1.5 in.

Weight: 2.7 oz

### HYFIN<sup>®</sup> TRAINER Twin Pack



- Unique adhesive gel formulation does not damage mannequins or leave excessive residue behind
- Packaged with two training chest seals in one unique perforated package allowing students to open only one seal at a time

#### SPECIFICATIONS

Packaged: L 7.5 in. x W 4.5 in. x D 0.25 in.

Chest Seal: H 6 in. x W 6 in. x D 0.1 in.

Weight: 2.5 oz

#### HYFIN TRAINER

##### ITEM #

10-0032

## WOUND MANAGEMENT SIMULATOR

- Rugged part-task trainer for learning essential life-saving skills of tourniquet application, wound packing and applying pressure bandages.
- Multiple wounding patterns to enhance training objectives
- Simulates severe laceration, penetrating trauma (entrance and exit wounds), along with the ability to apply tourniquets
- Proprietary silicon that is rugged and durable and will not break down due to environmental exposure
- Can be used in the classroom or in field training scenarios



*Portable training aid to help learn wound management techniques to control bleeding*

#### WOUND MANAGEMENT SIMULATOR

TONE	ITEM #
Light	93-0069
Dark	93-0410

#### SPECIFICATIONS

L 15.25 in. x Diameter 6 in.

Weight: 10.38 lb



## C-A-T® Combat Application Tourniquet Gen 7

**The Fastest, Safest, Most Effective Prehospital Field Tourniquet**

- **Single route buckle** decreases application time to help decrease blood loss, allows for slack removal for proper TQ application
- **Applied same way** to both upper and lower extremities
- True **one-handed tourniquet**
- **Scientifically proven** 100% effective in occluding blood flow in both upper and lower extremities by the USAISR and recommended by CoTCCC



online video available

scan to view>



### C-A-T TOURNIQUET

COLOR	STYLE	ITEM #
Black	Tactical	30-0001
Orange	Rescue	30-0023
Blue	Trainer	30-0033

### SPECIFICATIONS

Packaged: H 6.5 in. x W 2.4 in. x D 1.5 in.  
Device: L 37.5 in.  
Weight: 2.7 oz

## QuikClot® BLEEDING CONTROL DRESSING

- Advanced hemostatic agent combined with the simplicity of gauze available without medical device authorization (MDA)
- Utilizes the same inorganic component as QuikClot Combat Gauze® used by the US Military to control severe bleeding

### SPECIFICATIONS

**Z-fold Gauze** L 4 yd x W 3 in.  
Pkg: H 5 in. x W 5.5 in. x D 0.625 in.  
Weight: 0.8 oz

**Roll Gauze** L 4 ft x W 3 in.  
Pkg: H 5 in. x W 3.5 in. x D 0.5 in.  
Weight: 0.5 oz

### QUIKCLOT BCD

STYLE	ITEM #
Z-Fold	30-0161
Roll	30-0164



## NAR WOUND PACKING GAUZE

- Ultra compressed, soft, white 1-ply gauze
- Includes embedded X-Ray detectable line
- Z-fold packaging for controlled application as a complete unit or in a continuous linear feed
- Used as a packing gauze for wounds, as a backing gauze for hemostatic agents, with a compression bandage to control hemorrhage, or for minor wound bandaging

### SPECIFICATIONS

Gauze: W 4.5 in. x L 4.1 yd  
Packaged: L 3 in. x W 2 in. x D 1 in.  
Weight: 12 oz

### WOUND PACKING GAUZE

ITEM #  
30-0054

**Responder  
COMPRESSED GAUZE**

- Low-cube, ultra-compact sterile gauze
- Premium rolled 100% cotton gauze for maximum absorption

**SPECIFICATIONS**

Gauze: W 4.5 in. x L 4.1 yd  
 Packaged:  
 L 3 in. x W 2.5 in. x D 1 in.  
 Weight: 1.2 oz

**RESPONDER  
COMPRESSED GAUZE**

**ITEM #**

30-0182



**NAR S-ROLLED GAUZE**

- Unique inner dispenser, allows controlled linear feed application
- Premium 100% cotton

**SPECIFICATIONS**

Gauze: W 4.5 in. x L 4.1 yd  
 Packaged:  
 L 3 in. x W 2.5 in. x D 1 in.  
 Weight: 1.2 oz

**NAR  
S-ROLLED GAUZE**

**ITEM #**

30-0003

**RESPONDER  
ETD™ Emergency  
Trauma Dressing**

- Sterile pressure bandages that deliver extremely effective pressure and excellent absorption
- Suitable for applying immediate direct pressure
- Equipped with a non-adherent pad
- Easy-to-use hook & loop securing device with no confusing pressure bars or hooks
- Simple application with the Quick grip roll control system which prevents the dressing from unrolling if dropped.
- Rugged vacuum-sealed, low-cube packaging with easy-open tear notches



**RESPONDER EMERGENCY TRAUMA DRESSINGS**

STYLE	BANDAGE	PAD	WGT	PACKAGE	ITEM #
4 in. Flat	W 4 in. x L 60 in	4x6	2 oz	L 4.25 in. x W 2.75 in. x D 0.88 in.	30-0171
4 in. Flat (Mini)	W 4 in. x L 42 in.	4x6	2.4 oz	L 5 in. x W 3.5 in. x D 0.75 in.	30-0179
6 in. Rolled	W 6 in. x L 60 in.	6x6	2.9 oz	L 6 in. x W 2.5 in. x D 2.1 in.	30-0168

**RESPONDER  
Nitrile Gloves**



- Form-fitting, ambidextrous design with textured fingertips
- Thick palm ensures durability
- Textured for improved wet and dry gripping
- Longer and thicker than standard nitrile gloves

**SPECIFICATIONS**

Bag (25 pair): L 4.75 in. x W 6 in. x D 4 in.  
 Weight: 1 lb 1 oz


**RESPONDER NITRILE GLOVES**

SIZE	PKG   QTY	ITEM #
S	Bag   25 pair	ZZ-1059
M	Bag   25 pair	ZZ-1060
L	Bag   25 pair	ZZ-1043
XL	Bag   25 pair	ZZ-1061

\* also available in black at [NARescue.com](http://NARescue.com)

## HYFIN® VENT COMPACT Twin Pack

- Two Chest Seals packaged approximately 25% smaller than our standard HyFin®
- 3-channel pressure relief vents
- Two Chest Seals designed for the treatment of entry/exit or multiple penetrating wounds
- Large red tab designed for single step, peel and apply application and for burping the wound if necessary



### SPECIFICATIONS

Device (excludes tab):  
4.75 in. x 4.75 in.  
Wgt: 1.55 oz  
Packaged:  
H 3.88 in x W 6.6 in x D 0.13 in.

### HYFIN VENT COMPACT

ITEM #
10-0042

## NAR SURVIVAL BLANKET

Aluminized, reflective PET to help retain up to 90% of your body heat

### SURVIVAL BLANKET

ITEM #	NSN #
ZZ-0038	4240-01-697-8869

### SPECS

Unfolded: L 84 in. x W 52 in. | Weight: 17 oz



## Responder 7.25 in. TRAUMA SHEARS



Large, premium-grade Trauma Shears designed for safety and durability while cutting through tough materials

### SPECIFICATIONS

L 7.25 in. x W 3.7 in. x D 0.4 in.  
Weight: 2 oz

### SHEARS

ITEM #
ZZ-1019

## Responder QUIKLITTER

- Multi-use, low cost disposable evacuation platform for effective emergency casualty transport & transfer
- Constructed of a latex free water-resistant, woven polyethylene with signature red nylon webbing & stitching
- Ten (10) carry handles with four (4) on each side and one (1) larger handle at either end of the litter to provide multiple carry points.



### QUIKLITTER

ITEM #
60-0080

### SPECIFICATIONS

Payload 500 lb.  
Litter Bed: L 74 in. x W 28 in.  
Weight: 13.3 oz

## Emergency SURVIVAL WRAP

- Lightweight, durable, laminated polyethylene
- Reversible Silver/Orange or OD Green



### SURVIVAL WRAP

COLOR	ITEM #
Orange/Silver	ZZ-0834
OD Green/Silver	ZZ-1390

### SPECIFICATIONS

Packaged:  
L 5 in. x W 3.75 in.  
Open Wrap:  
L 96 in. x W 60 in.  
Weight: 3.7 oz

PUBLIC ACCESS  
**BLEEDING CONTROL**  
INDIVIDUAL KITS | PACKS | STATIONS



COMMUNITY  
PREPAREDNESS



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277